

Cells with or without CGK733 treatment	Cells with abnormal karyotype/ total number of cells analyzed	Karyotype
Near replicative senescent cells	0/50	46, XY
Replicative senescent cells + CGK733	0/50	46, XY
BJ cells infected with vector	0/50	46, XY
TRF2 ^{ΔBΔM} -induced senescent cells + CGK733	3/50	46, XY [47] ^a 46, XY, t(7;14) [1] 47, XYY [1] 47, XY, +10 [1]

^aBrackets [], number of cells.

Supplementary Figure 3. Karyotypes of the senescent cells induced to proliferate by CGK733. For cells treated with CGK733, karyotyping analysis was performed at 2 weeks after CGK733 treatment.