


Sample Information

Raw File Name	D:\Data\4713\4713_27.raw
Instrument Method	C:\Xcalibur\methods\UltiMate\NoFAIMS_Intact_Protein\Direct_Injection_MS1_IT_7K_RF60_35min.meth
Vial	RC3
Injection Volume (µL)	1
Sample Weight	0
Sample Volume (µL)	0
ISTD Amount	0
Dil Factor	1

Chromatogram Parameters


Use Restricted Time	True
Time Limits	15.000 - 24.984 minutes
Scan Range	558 - 930
m/z Range	600 - 2000
Chromatogram Trace Type	TIC
Sensitivity	High
Rel. Intensity Threshold (%)	5


Chromatogram



Main Parameters (ReSpect™)	
Deconvolution Results Filter	
Output Mass Range	22500 - 35000
Deconvoluted Spectra Display Mode	Isotopic Profile (new)
Charge State Distribution	
Deconvolution Mass Tolerance	30 ppm
Choice of Peak Model	
Choice of Peak Model	Intact Protein
Resolution at 400 m/z	
Raw File Specific	2000
Generate XIC for Each Component	
Calculate XIC	True
Advanced Parameters (ReSpect™)	
Charge State Distribution	
Model Mass Range	8000 - 70000
Charge State Range	7 - 100
Minimum Adjacent Charges (low & high model mass)	4 - 4
Noise Parameters	
Rel. Abundance Threshold (%)	0
Deconvolution Quality	
Quality Score Threshold	0
Choice of Peak Model	
Target Mass	28000 Da
Peak Model Parameters	
Number of Peak Models	1
Left/Right Peak Shape	2:2
Peak Filter Parameters	
Peak Detection Minimum Significance Measure	1 Standard Deviations
Peak Detection Quality Measure	95%
Specialized Parameters	
Peak Model Width Factor	1
Intensity Threshold Scale	0.01
Deconvolution Parameters	
Noise Compensation	True
Charge Carrier	H
Negative Charge	False
Source Spectra Parameters	
Source Spectra Method	Average Over Selected Retention Time
RT Range	16.000 - 22.000 minutes

4713_27 #595-818 RT:16.000-22.000 AV:224
F:ITMS + p NSI Full ms [600.0000-2000.0000]





ReSpect Masses Table													
Row Number	Average Mass	Intensity	Relative Abundance	Fractional Abundance	Score	Number of Charge States	Charge State Distribution	Mass Std Dev	PPM Std Dev	Delta Mass	Start Time (min)	Stop Time (min)	Apex RT
1	25416.26	815090.63	100.00	23.06	65.78	15	15 - 29	1.67	65.53	0.00	16.000	22.000	17.890
2	25372.77	767286.69	94.14	21.71	38.13	8	28 - 35	1.11	43.66	-43.49	16.000	22.000	18.610
3	25456.06	636106.81	78.04	17.99	39.66	8	22 - 29	1.76	69.21	39.80	16.000	22.000	18.590
4	24767.06	182746.95	22.42	5.17	20.61	4	27 - 30	0.96	38.79	-649.20	16.000	22.000	18.610
5	24437.86	173108.19	21.24	4.90	27.50	6	22 - 27	1.30	53.28	-978.40	16.000	22.000	18.160
6	26958.56	164699.78	20.21	4.66	18.17	4	32 - 35	2.51	93.13	1542.30	16.000	22.000	18.610
7	22993.64	141204.67	17.32	3.99	18.98	4	28 - 31	1.29	56.20	-2422.62	16.000	22.000	18.240
8	27193.89	107299.01	13.16	3.04	17.75	4	23 - 26	1.85	68.04	1777.63	16.000	22.000	18.830
9	25983.92	101641.03	12.47	2.88	16.69	4	22 - 25	1.95	74.91	567.66	16.000	22.000	18.590
10	28998.29	93763.45	11.50	2.65	19.71	4	26 - 29	0.36	12.27	3582.03	16.000	22.000	18.800
11	26065.04	77241.05	9.48	2.19	26.01	5	20 - 24	2.43	93.08	648.78	16.000	22.000	18.610
12	28374.89	70934.79	8.70	2.01	18.46	4	22 - 25	2.49	87.59	2958.63	16.000	22.000	18.800
13	25371.09	61942.67	7.60	1.75	19.07	4	17 - 20	1.77	69.64	-45.17	16.000	22.000	18.860
14	32635.55	41435.00	5.08	1.17	15.70	4	29 - 32	3.50	107.12	7219.29	16.000	22.000	18.690
15	29325.76	28532.39	3.50	0.81	14.93	4	30 - 33	2.49	84.89	3909.50	16.000	22.000	18.100
16	26105.05	25034.94	3.07	0.71	18.97	4	14 - 17	1.90	72.60	688.79	16.000	22.000	18.750
17	31886.56	24281.38	2.98	0.69	19.91	4	16 - 19	3.06	96.00	6470.30	16.000	22.000	17.590
18	26023.38	22686.74	2.78	0.64	24.76	5	14 - 18	2.08	79.78	607.12	16.000	22.000	17.830