


Sample Information

Raw File Name	D:\Data\4713\4713_08.raw
Instrument Method	C:\Xcalibur\methods\UltiMate\NoFAIMS_Intact_Protein\Direct_Injection_MS1_IT_7K_RF60_35min.meth
Vial	RA8
Injection Volume (µL)	1
Sample Weight	0
Sample Volume (µL)	0
ISTD Amount	0
Dil Factor	1

Chromatogram Parameters


Use Restricted Time	True
Time Limits	15.000 - 24.984 minutes
Scan Range	558 - 930
m/z Range	600 - 2000
Chromatogram Trace Type	TIC
Sensitivity	High
Rel. Intensity Threshold (%)	5


Chromatogram



Main Parameters (ReSpect™)	
Deconvolution Results Filter	
Output Mass Range	22500 - 35000
Deconvoluted Spectra Display Mode	Isotopic Profile (new)
Charge State Distribution	
Deconvolution Mass Tolerance	30 ppm
Choice of Peak Model	
Choice of Peak Model	Intact Protein
Resolution at 400 m/z	
Raw File Specific	2000
Generate XIC for Each Component	
Calculate XIC	True
Advanced Parameters (ReSpect™)	
Charge State Distribution	
Model Mass Range	8000 - 70000
Charge State Range	7 - 100
Minimum Adjacent Charges (low & high model mass)	4 - 4
Noise Parameters	
Rel. Abundance Threshold (%)	0
Deconvolution Quality	
Quality Score Threshold	0
Choice of Peak Model	
Target Mass	28000 Da
Peak Model Parameters	
Number of Peak Models	1
Left/Right Peak Shape	2:2
Peak Filter Parameters	
Peak Detection Minimum Significance Measure	1 Standard Deviations
Peak Detection Quality Measure	95%
Specialized Parameters	
Peak Model Width Factor	1
Intensity Threshold Scale	0.01
Deconvolution Parameters	
Noise Compensation	True
Charge Carrier	H
Negative Charge	False
Source Spectra Parameters	
Source Spectra Method	Average Over Selected Retention Time
RT Range	16.000 - 22.000 minutes

4713_08 #596-819 RT:16.000-22.000 AV:224
F:ITMS + p NSI Full ms [600.0000-2000.0000]





ReSpect Masses Table													
Row Number	Average Mass	Intensity	Relative Abundance	Fractional Abundance	Score	Number of Charge States	Charge State Distribution	Mass Std Dev	PPM Std Dev	Delta Mass	Start Time (min)	Stop Time (min)	Apex RT
1	24185.57	7876414.50	100.00	21.47	34.88	7	19 - 25	1.26	52.08	0.00	16.000	22.000	17.720
2	27010.35	5086966.50	64.58	13.87	23.54	4	28 - 31	1.00	36.92	2824.78	16.000	22.000	17.460
3	22941.15	3349989.50	42.53	9.13	16.17	4	21 - 24	1.59	69.40	-1244.42	16.000	22.000	17.240
4	22612.91	3270978.00	41.53	8.92	21.42	4	24 - 27	1.92	85.11	-1572.67	16.000	22.000	17.130
5	22624.09	2026613.13	25.73	5.52	15.21	4	29 - 32	2.26	99.78	-1561.48	16.000	22.000	17.030
6	32248.49	1972378.63	25.04	5.38	21.16	4	31 - 34	1.60	49.62	8062.92	16.000	22.000	17.480
7	33050.26	1943441.50	24.67	5.30	19.50	4	31 - 34	2.46	74.49	8864.69	16.000	22.000	17.210
8	29511.66	1646949.88	20.91	4.49	16.21	4	29 - 32	1.98	67.00	5326.08	16.000	22.000	17.350
9	26853.65	1286400.25	16.33	3.51	27.13	5	18 - 22	1.70	63.13	2668.08	16.000	22.000	17.370
10	31667.22	1261647.25	16.02	3.44	25.31	5	20 - 24	1.47	46.37	7481.65	16.000	22.000	17.640
11	26354.16	1141496.63	14.49	3.11	27.59	6	19 - 24	0.90	34.07	2168.58	16.000	22.000	17.480
12	27500.33	872926.38	11.08	2.38	20.39	4	21 - 24	2.91	105.74	3314.76	16.000	22.000	17.130
13	27972.02	749320.50	9.51	2.04	19.60	4	23 - 26	2.20	78.64	3786.45	16.000	22.000	17.480
14	27737.08	666703.81	8.46	1.82	20.87	4	22 - 25	1.63	58.84	3551.51	16.000	22.000	17.240
15	26087.63	528168.75	6.71	1.44	20.76	4	17 - 20	2.25	86.26	1902.06	16.000	22.000	17.400
16	27008.90	516173.53	6.55	1.41	16.73	4	20 - 23	1.94	71.86	2823.33	16.000	22.000	17.210
17	27317.74	514621.06	6.53	1.40	20.77	4	19 - 22	1.55	56.80	3132.17	16.000	22.000	17.110
18	30119.15	500826.63	6.36	1.37	21.64	5	16 - 20	1.90	63.10	5933.58	16.000	22.000	17.560
19	30691.04	408929.34	5.19	1.11	23.38	5	22 - 26	1.57	51.02	6505.47	16.000	22.000	17.210
20	28986.14	390871.66	4.96	1.07	20.70	4	19 - 22	2.46	85.03	4800.57	16.000	22.000	17.190
21	27991.22	282018.09	3.58	0.77	20.09	4	18 - 21	2.18	77.97	3805.65	16.000	22.000	17.210
22	25541.28	262586.34	3.33	0.72	17.51	4	14 - 17	1.83	71.52	1355.71	16.000	22.000	17.370
23	23923.74	124868.29	1.59	0.34	21.81	5	16 - 20	1.64	68.66	-261.83	16.000	22.000	17.480