
Supplementary information


Lesion recognition by XPC, TFIIH and XPA in DNA excision repair

In the format provided by the
authors and unedited


Supplementary Table 1. Oligonucleotide sequences used in EMSA, cryoEM and dual incision assay

Oligo Name	Sequence (5'→3')
Top56_iCy5	5'-AGGTAGTCACAGCTGATTGCGCTGAGGAT- iCy5 -TAGACGTGCGATATCTATGCCACTT
Top56_2AP	5'-AGGTAGTCACAGCTGATTGCGCTGAGGAT- ApAp -TAGACGTGCGATATCTATGCCACTT
Bot56	5'-AAGTGGCATAGATATCGCACGTCTA TT TATCCTCAGCGCAATCAGCTGTGACTACCT
Top_F1_57	5'-ACTGGGAACAGCTGCGGTCCGTTGACGACATACAGCTCAACGATGCCACTTCAGCAT
Top_F1_47	5'-GCTGCGGTCCGTTGACGACATACAGCTCAACGATGCCACTTCAGCAT
Top_F1_37	5'-GTTGACGACATACAGCTCAACGATGCCACTTCAGCAT
Top_F2_iCy5	/5Phos/TGCACTCATGTTCTG- iCy5 -CACTAGTACTACGACTACGC
Top_F3_10_5P	/5Phos/CTCGGCTCGA
Top_F3_20_5P	/5Phos/CTCGGCTCGAGTGTTCGATC
Top_F1_53	5'-GCTGCGGTCCGTTGACGACATACAGCTCAACGATGCCACTTCAGCATTGCACT
Top_F2_AP	³² p-5'-CATGTTCTG- ApAp -CACTAGTAC
Top_F3_21_5P	/5Phos/TACGACTACGCCTCGGCTCGA
Cy5_72/20_top	Top_F1_57+Top_F2_iCy5
Cy5_62/30_top	Top_F1_47+Top_F2_iCy5+Top_F3_10_5P
Cy5_52/40_top	Top_F1_37+Top_F2_iCy5+Top_F3_20_5P
AP_62/30_top	Top_F1_53+Top_F2_AP+Top_F3_21_5P
Bot_F1_74	5'-CGAACACTCGAGCCGAGGCGTAGTCGTAGTACTAGTGTTTCAGAACATGAGTGCAATGCTGA AGTGGCATCGTTG (a template to ligate top_F1, F2 and F3)
mm94_Bot_72/20	5'-GCGTAGTCGTAGTACTAGTG TT TCAGAAATATGAGTACAATACTGAAGTGGCATCGTTGAGCTG TATGTCGTCAACGGACCGCAGCTGTTCCCACT
mm94_Bot_62/30	5'-TCGAGCCGAGGCGTAGTCGTAGTACTAGTG TT TCAGAAATATGAGTACAATACTGAAGTGGCAT CGTTGAGCTGTATGTCGTCAACGGACCGCAGC
mm94_Bot_52/40	5'-GATCGAACACTCGAGCCGAGGCGTAGTCGTAGTACTAGTG TT TCAGAAATATGAGTACAATACT GAAGTGGCATCGTTGAGCTGTATGTCGTCAAC
Annealed Substrate (Diagram)	
Cy5-DNA	Top56_iCy5 5' _____ Cy5 _____
	Bot56 3' _____ 29 _____ TT _____ 25 _____
AP-DNA	Top56_2AP 5' _____ AP _____
	Bot56 3' _____ 29 _____ TT _____ 25 _____
Cy5- or AP-DNA in dual incision	Top94 (F1+F2 (+F3)) 5' _____ XX _____
	mm94_Bot (3 variation) 3' _____ 52 / 62 / 72 _____ TT _____ 40 / 30 / 20 _____


Supplementary Figure 1 Uncropped SDS gels related to Fig. 1b




Supplementary Figure 2 Uncropped EMSA gels related to Fig. 1c



Supplementary Figure 3 Uncropped EMSA gels related to ED Fig. 3a



Supplementary Figure 4 Gels of DNA dual incision assay used to generate Fig. 4c



NER prod (%)							
	Cy5-DNA			AP-DNA			
min	run 1	run 2	run 3	run 1	run 2	run 3	
0	0	0	0	0	0	0	
5	14.48	14.44	12.21	1.79	3.63	1.18	
10	27.41	24.72	27.72	4.42	7.94	2.46	
20	41	45.15	47.08	9.27	14.92	7	
30	50.77	60.46	51.03	14.8	20.28	9.38	
60	63.16	75.52	62.46	24.47	28.53	13.6	

Supplementary Figure 5 Gels of DNA dual incision assay used to generate ED Fig. 10d

