
Supplementary information

Quality control of glycogen through direct ubiquitylation by RNF213

In the format provided by the
authors and unedited

Fig. 2g




Fig. 3d




Fig. 3f







Fig. 2i



Extended Data Fig. 1a




Extended Data Fig. 1h





note: same as GYS1 (above), but rotated to ensure bands of interest are straight

TMR



Extended Data Fig. 1i





Extended Data Fig. 2a,b




Extended Data Fig. 3j




Extended Data Fig. 3l




Extended Data Fig. 4c




Extended Data Fig. 4a




Extended Data Fig. 4b




Extended Data Fig. 7b



Extended Data Fig. 8a



Extended Data Fig. 8h



Extended Data Fig. 9f,g

(note: same blot, just rotated to ensure bands of interest are straight)

