




Note on uncropped blots: In our experimental protocol, the PVDF membranes were cut horizontally into strips prior to blocking and antibody incubation based on the pre-stained molecular weight markers. This standard procedure was performed to allow simultaneous probing of multiple target proteins with different molecular weights from a single gel and to conserve primary antibodies. Therefore, full-length (whole-gel) blots are not available. The images provided below represent the original, entirely uncropped, and unprocessed scans of the individual cut membrane strips.





Eif4e2 (28kDa)

← 37 kDa
← 30 kDa
← 28 kDa
← 16 kDa

