

MSK Skin Tone Protocol Checklist:

- ☐ All assessments were performed in a single room with standardized environment (humidity/temperature) and consistent ambient lighting. Patients acclimated to the environment for 5 minutes before any color assessments.
- ☐ Demographic information
 - ☐ Race (White, Black, Asian, American Indian/Alaskan, Hawaiian/Pacific, Other)
 - ☐ Ethnicity (Hispanic vs Non-Hispanic)
- ☐ Fitzpatrick questionnaire
 - ☐ Tanning propensity (Do not tan, tan poorly, tan after burn, tan easily, tans easily darkly, always tan darkly)
 - ☐ Burning propensity (Always burn, easily burn, burn then tan, burn minimally, rarely burn, never burn)
 - ☐ Natural Hair Color
 - ☐ Natural Eye Color
- ☐ Patient Self perceived Pantone and Monk at inner wrist
- ☐ Up to 10 lesions and non-lesion anatomic site were marked with a sticker
 - ☐ Constitutive sites: inner upper arm, ventral forearm, abdomen, lower back
 - ☐ Facultative sites: forehead, dorsal forearm, chest, upper back, shin, calf, sole
- ☐ Total Body Photography in white light and cross polarized light
- ☐ Pantone perception by 2 raters, Monk perception by 2 raters
- ☐ Pantone color card iPhone App at inner wrist (x2), outer forearm, inner upper arm
- ☐ Dermoscopic images with camera #1 of each lesion & anatomic site in non-contact white and cross polarized and contact white and cross polarized
- ☐ Dermoscopic images with camera #2 at lesion #1, outer forearm, inner wrist
- ☐ Colorimeter readings at inner wrist (x3) forearm (x3), upper inner arm(x3), skin around lesion 1 (x3)



