

PCR primers used for preparation of FISH probes for reconstructed tandem repeats

Repeat	Primers (5'->3')
TR-1	TTCAACAACGGCAAAACATC, CTATATACTTTTAAGATCAATGG
TR-2	GAAGGCCATCCAAGTTTCAT, ATGACCTATAATCTTTCACC
TR-3	ATCTCATTTATATAATGAGG, AATATACATATKACAACCTC
TR-4	GTCTTAGATTCAATTAAGGT, AAAATTGAATGAGACTCTTC
TR-5	CGAGTGTTGAAAATCAATAC, TCATTTCTTTGTCAATACTC
TR-6	TGATTTCAATGAATAACTCG, ATATGATGTATGGGTATTC
TR-7	GTATTTGTATACATGTTTTACC, AGGATAAATGTGGTACCTAC
TR-8	TTGCCAGCGGAATTATCAC, AATAACACCAAACACCACTC
TR-9	ATACCATAATAYAATGTGAATAT, CAAATTTTAGGTTTCTTCATC
TR-10	TTATTATCAGAACCATGAAC, TGTAAGAAAGACTTCTTTATGG
TR-11	TRTTACTTAGCTACTTAATG, CCATCCTATAATGTTCAAAG
TR-12	AATAGCCTTTCCTTAGTGTC, GATTTAATCGACATTGGAAG
TR-14	CAACCCAAAACAAAGAAATTAC, GGTTTTGGCTCTTCGTAATC
TR-17	AGTCAATTACGGTGTAGGTG, TACATTGATTTGTGTGCAAG