

Table 1. Number of CFU recovered on BAP after co-culturing infected cells with fresh host cells

Addition of new host cells		Number of CFU recovered on BAP after 24 hour co-culture		
		iEC	iSMC	iKB
Co-culture	Medium only	0	0	33 +/- 14
	KB	72 +/- 8	400 +/- 12	2169 +/- 107
	SMC	26 +/- 6	210 +/- 18	191 +/- 3
Transwell System	Medium only	0	0	0
	KB	7 +/- 3	64 +/- 11	361 +/- 19
	SMC	2 +/- 1	30 +/- 5	79 +/- 15