

A**IVD**

4 °C 30 °C

Commercial name, Malaria Antigen Pf/Pan (HRP2/ pLDH) RDT**REF**

XXXXXXXXXX



xx

Rapid test for the antigen detection of malaria (*Plasmodium x*)**Test rapide de détection d'antigène du paludisme (*Plasmodium x*)****Prueba rápida para detección de antígeno de malaria (*Plasmodium x*)****Teste rápido para detecção de antígeno da malária (*Plasmodium x*)**

		LOT	XXXXXXXX
			YYYY-MM
	Name, address, website, telephone and/or fax number	CE	EC REP Name and address

B**Commercial name, Malaria Antigen Pf/Pan (HRP2/ pLDH) RDT****REF**

XXXXXXXXXX

Content :	Cassettes Alcohol swabs Lancets Specimen transfer devices Buffer bottle Instructions for use	 LOT	xx XXXXXXXX
Required but not provided	Timer Gloves Biosafety sharps container Waste container		YYYY-MM
		4 °C 30 °C	