


A



B

	Forward primer	Reverse primer
(a) clone4-35	5'-CCCAACCACTCCCCAAACAGGCTGG-3'	5'-TTACAGCCGTTGCTCTCCTT-3'
(b) clone41	5'-GCCATGCTGATGTCAGAGAA-3'	5'-GATCCCCAGTCAGGATGAGA-3'
(c) 143-145	5'-TTGGTCCTGGGTGCTCAAAT-3'	5'-TTTCCCAAGAGTACGGCAGT-3'
(d) iPCR145	5'-GAATCCCCATCTTAGCATCTAAGGGATT-3'	5'-CTGGAAATACTGTTCTTGAGGTCATGGT-3'
(e) AKF1-10	5'-GGCCATAGACCCCTCTTCTCAGTAA-3'	5'-CTCGTGCATTTGGCTGGGGAGTGATTC-3'
(f) AK126481	5'-GCTCAGTGCCAGCTGCTTAAAAAT-3'	5'-TAGTCCTGGGCAGACCAGTTTCTATC-3'
LOC3	5'-GACCAGACCCCAGGAAAGA-3'	5'-TTACAGCCGTTGCTCTCCTT-3'
LOC10	5'-AGCAAGAACTCTGGAGAAGCA-3'	5'-GAGAGGCGTGGGTGAGAG-3'
GAPDH(for real-time)	5'-CCACATCGCTCAGACACCAT-3'	5'-GCAACAATATCCACTTTACCAGAGTTAA-3'
GAPDH(for cloning)	5'-CCACCCATGGCAAATTCCATGGCA	5'-TCTAGACGGCAGGTCAGGTCCACC-3'