

Precisions and repeatability

Peak No.	Saponins	Intra-day Precision (%)	Inter-day Precision (%)	Repeatability
1	Notoginsenoside R ₁	2.2	3.4	2.4
2	Ginsenoside Rg ₁	1.0	1.4	0.6
3	Ginsenoside Re	7.7	5.1	0.88
4	Notoginsenoside R ₄	2.2	4.1	3.6
5	Ginsenoside Rf	2.7	0.5	5.4
6	Notoginsenoside Fa	1.9	5.9	2.5
7	Notoginsenoside I	2.3	3.6	2.3
8	SC1	0.9	3.3	3.4
9	Ginsenoside Rb ₁	0.5	2.9	1.1
10	Notoginsenoside Fc	2.7	3.3	5.1
11	Ginsenoside Rg ₂	1.3	2.8	3.3
12	Ginsenoside Rh ₁	0.7	2.8	2.1
13	Ginsenoside Rb ₂	0.9	1.6	5.4
14	Ginsenoside F ₁	3.3	2.6	7.9
15	Ginsenoside Rd	1.4	3.8	0.7
16	Notoginsenoside K	1.0	1.0	2.9
17	Notoginsenoside T ₅ /Unkown	3.7	0.8	9.1
18	Unkown	0.7	4.6	7.3
19	Notoginsenoside T ₅ /Unkown	3.5	5.5	2.4
20	Unkown	1.8	2.9	7.7
21	Ginsenoside Rk ₃	2.2	3.3	2.4
22	Ginsenoside Rh ₄	0.7	3.6	2.2
23	20(S)-Rg ₃	2.1	2.2	3.0
24	20(R)-Rg ₃	1.9	2.4	4.7
25	Ginsenoside F ₂	6.3	8.6	3.8
26	Ginsenoside Rk ₁	2.6	3.0	4.5
27	Ginsenoside Rg ₅	5.9	3.4	4.4