

Recovery

Analytes	Added (µg/ml)	found(µg/ml, n=6)	Recovery (% , n=6)	RSD (% , n=6)
Notoginsenoside R ₁	237.50	239.96	101.0	4.7
Ginsenoside Rg ₁	658.75	660.75	100.3	7.3
Ginsenoside Re	302.50	325.44	107.6	3.2
Ginsenoside Rb ₁	506.25	505.05	99.8	4.5
Ginsenoside Rg ₂	25.00	23.92	95.7	4.0
Ginsenoside Rh ₁	33.75	31.97	94.7	1.4
Ginsenoside Rb ₂	31.25	28.94	92.6	2.6
Ginsenoside Rd	283.75	296.23	104.4	3.2
Ginsenoside 20(S)-Rg ₃	23.63	24.76	104.8	3.4
Ginsenoside 20(R)-Rg ₃	7.00	7.08	101.1	9.1
Ginsenoside 20(R)-Rg ₃ [*]	7.00	6.93	99.0	7.4

^{*} The recovery and RSD results were obtained by using the calibration curve for 20(S)-Rg₃ as the substitutive one for 20(R)-Rg₃.