





45 Dentin Specimens

Half of the 45 Dentin specimens were covered with a self-etching adhesive system
R: Reference site

The specimens were placed in the intraoral appliances

Control group; (15 specimens)
- 10 ml distilled water for 60 sec. daily for 3 weeks

Wash-out one week

250 ppm Fluoride group; (15 specimens)
- 10 ml mouth rinse for 60 sec. daily for 3 weeks

Wash-out one week

5000 ppm Fluoride group; (15 specimens)
- 1:3 slurry 10 ml, 60 sec. daily for 3 weeks

TMR Analysis