

Table S2. Homeolog-specific primers used for real-time qPCR analysis.

Homeolog	Forward (5'→3')	Reverse (5'→3')	Amplicon (bp)
<i>PSY-A1</i>	GAATGAAACATGGCAGTGTG	AGAACTTCATGTCATGCCTA	120
<i>PSY-B1</i>	GATGAAGATGGTCAGGTTAG	CAGAACTTCATGTTATGCATG	161
<i>PSY-D1</i>	AAAGAAACAGGGCGGTGTC	AGCCCTTCGTGTCATGGATT	111
<i>LCYe-A</i>	GTTTCGACGCTTTCGTCGGT	TCTCGCAGGCTAGTCAAATG	175
<i>LCYe-B</i>	GGATATCACCTGCCATGTG	CACGCTCCCACTTTTGGCA	144
<i>LCYe-D</i>	GGATATCACCTGCCATGTG	AACCTAAAATTAAACCATGGCA	102
<i>CCD-A1</i>	TGTCAACGAGGAACAACCTGG	GCAATGCTACTAATTTAACGC	85
<i>CCD-B1</i>	AGAAGCAGCTTGAAGTCGGT	GATAGCCGTCATCTTCTTGA	126
<i>CCD-D1</i>	TGTCAACGAGGAACAACCTGG	TCTGCAGCATGCAGACCAC	98
<i>CCD-A4</i>	CCGAGCTCCGGTCACAGC	ATACATCGCGCGATCAATCC	89
<i>CCD-B4</i>	TCCACGGCCTATTCGTCACT	CTATTATACATCGTGGCGCG	121
<i>CCD-D4</i>	CCGAGCTCCGGTCACAGC	TCTACTATTGTACATTGTCGCA	103