







UV-vis analysis for the detection of 500 PFU of RSV at 30, 60 and 120 min. using functionalized copper (A, B and C) silver (D, E and F) and gold nanoparticles (G, H and I) with 3 replicates.




UV-vis analysis for the detection of 500 PFU of RSV at 30, 60 and 120 min. using copper (A, B and C) silver (D, E and F) and gold nanoparticles (G, H and I) with 3 repetitions.




UV-vis analysis for the detection of 1000 PFU of RSV at 30, 60 and 120 min. using functionalized copper (A, B and C) silver (D, E and F) and gold nanoparticles (G, H and I) with 3 repetitions.



UV-vis analysis for the detection of 1000 PFU of RSV at 30, 60 and 120 min using copper (A, B and C) silver (D, E and F) and gold nanoparticles (G, H and I) with 3 repetitions.



UV-vis analysis for the detection of 2000 PFU of RSV at 30, 60 and 120 min. using functionalized copper (A, B and C) silver (D, E and F) and gold nanoparticles (G, H and I) with 3 repetitions.



UV-vis analysis for the detection of 2000 PFU of RSV at 30, 60 and 120 min. using copper (A, B and C) silver (D, E and F) and gold nanoparticles (G, H and I) with 3 repetitions.