

Gene	Direction	Primer sequence (5' to 3')
TNF- α	Forward	5'- CCTGTAGCCCACGTCGTAG-3'
	Reverse	5'- GGGAGTAGACAAGGTACAACCC -3'
IL-1 β	Forward	5'-GAAATGCCACCTTTTGACAGTG -3'
	Reverse	5'-TGGATGCTCTCATCAGGACAG-3'
iNOS	Forward	5'-GGAGTGACGGCAAACATGACT-3'
	Reverse	5'-TCGATGCACAACTGGGTGAAC-3'
TGF- β	Forward	5' -TCTGCATTGCACTTATGCTGA-3'
	Reverse	5' -AAAGGGCGATCTAGTGATGGA-3'
Arg1	Forward	5' -CAAGACAGGGCTCCTTTCAG-3'
	Reverse	5' -TGGCTTATGGTTACCCTCCC-3'
IL-10	Forward	5' -AGGCGCTGTCATCGATTCT-3'
	Reverse	5' -TGGAGTCCAGCAGACTCAAT-3'
GAPDH	Forward	5' -TGGCCTTCCGTGTTCTTAC-3'
	Reverse	5' - GAGTTGCTGTTGAAGTCGCA-3'

Table S1. Primers used for quantitative real time PCR assays.