

Mitofusin 1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 1/ set A
cal. mean

10981867,5

cal. mean

25207886

	calibrator		11430450	1,04	26237938	1,04	1,00
A1	H1/B	Control 17/B	9591795	0,87	25525559	1,01	0,86
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	20146728	0,80	#####
A4	Sample X/T	Sample X/T	excl.	#####	19557452	0,78	#####
A5	H3/B	Control 3/B	9857655	0,90	20382535	0,81	1,11
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	10237275	0,93	20062467	0,80	1,17
A8	CR4/B/D7	ICU 4/B/D7	8426700	0,77	19913007	0,79	0,97
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	9640485	0,88	16150422	0,64	1,37
A12	CR8/B/D7	ICU 11/B/D7	10728945	0,98	17634541	0,70	1,40
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	9848925	0,90	17586742	0,70	1,29
A16	CR10/B/D7	ICU 13/B/D7	9398565	0,86	18543098	0,74	1,16
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	9021330	0,82	20294868	0,81	1,02
A20	CR11/B/D7	ICU 9/B/D7	10016415	0,91	21152021	0,84	1,09
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		10533285	0,96	24177834	0,96	1,00

Mitofusin 1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

10873196

cal. mean

28899095,5

	calibrator		10859124	1,00	27854350	0,96	1,04
B1	H4/B	Control 12/B	7693488	0,71	28818080	1,00	0,71
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	7412370	0,68	29889932	1,03	0,66
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	6837078	0,63	23959672	0,83	0,76
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	29444064	1,02	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	27218925	0,94	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25503469	0,88	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22796640	0,79	#####
B11	CR9/B/D0	ICU 3/B/D0	7719515	0,71	26684496	0,92	0,77
B12	CR9/B/D7	ICU 3/B/D7	6568250	0,60	27844674	0,96	0,63
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	4411132	0,41	29392896	1,02	0,40
B16	CR12/B/D7	ICU 20/B/D7	6324978	0,58	29492584	1,02	0,57
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	8398324	0,77	28194210	0,98	0,79
B20	CR13/B/D7	ICU 22/B/D7	9962825	0,92	25151087	0,87	1,05
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		10887268	1,00	29943841	1,04	0,97

Mitofusin 1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C
cal. mean
10128671
cal. mean
34079936

	calibrator		10666480	1,05	34138962	1,00	1,05
C1	H15/B	Control 2/B	7824968	0,77	29730870	0,87	0,89
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	8223374	0,81	30098518	0,88	0,92
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	9694316	0,96	20994664	0,62	1,55
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	3264620	0,32	22660589	0,66	0,48
C8	CR3/B/D7	ICU 5/B/D7	3680046	0,36	18025875	0,53	0,69
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	2987240	0,29	19994035	0,59	0,50
C12	CR5/B/D7	ICU 8/B/D7	2749512	0,27	17391984	0,51	0,53
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	5866150	0,58	22085112	0,65	0,89
C16	CR7/B/D7	ICU 6/B/D7	6670874	0,66	18388750	0,54	1,22
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	6890156	0,68	24852613	0,73	0,93
C20	CR20/B/D7	ICU 7/B/D7	8036200	0,79	25778115	0,76	1,05
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		9590862	0,95	34020910	1,00	0,95

Mitofusin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	9740400	cal. mean	35190284		
calibrator			10686850	1,10	37167787	1,06	1,04
D1	H2/B	Control 4/B	10420656	1,07	30949908	0,88	1,22
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	8056348	0,83	31992571	0,91	0,91
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	12798994	1,31	24406066	0,69	1,89
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	7979298	0,82	22372376	0,64	1,29
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	10037660	1,03	18229138	0,52	1,99
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	8955464	0,92	17679473	0,50	1,83
D12	CR2/B/D7	ICU 14/B/D7	8405764	0,86	18461459	0,52	1,64
D13	CR2/B/D180	ICU 14/B/D180	8926760	0,92	19948365	0,57	1,62
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	6749442	0,69	22814050	0,65	1,07
D18	CR17/B/D7	ICU 23/B/D7	4939020	0,51	24909295	0,71	0,72
D19	CR17/B/D180	ICU 23/B/D180	6052726	0,62	26712850	0,76	0,82
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			8793950	0,90	33212781	0,94	0,96

Mitofusin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	14986089	cal. mean	13612646,5		
calibrator			15077844	1,01	13056250	0,96	1,05
E1	H12/B	Control 7/B	11489400	0,77	15730101	1,16	0,66
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	5244390	0,35	17836110	1,31	0,27
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	10341720	0,69	17137798	1,26	0,55
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	15340140	1,02	16121115	1,18	0,86
E8	CR6/B/D7	ICU 21/B/D7	12693420	0,85	12428020	0,91	0,93
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	13375755	0,89	12644640	0,93	0,96
E12	CR14/B/D7	ICU 10/B/D7	12600135	0,84	13423200	0,99	0,85
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####
E19	CR19/B/D0	ICU 15/B/D0	11212740	0,75	18376065	1,35	0,55
E20	CR19/B/D7	ICU 15/B/D7	9039825	0,60	18750060	1,38	0,44
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			14894334	0,99	14169043	1,04	0,95

Mitofusin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	17262949	cal. mean	24950732		
calibrator			15682642	0,91	25509639	1,02	0,89
F1	H9/B	Control 9/B	13160462	0,76	32156010	1,29	0,59
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	12082498	0,70	31970444	1,28	0,55
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	13886894	0,80	25759862	1,03	0,78
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	12621756	0,73	24703341	0,99	0,74
F8	CR1/B/D7	ICU 1/B/D7	16108234	0,93	22559436	0,90	1,03
F9	CR1/B/D180	ICU 1/B/D180	11964094	0,69	26013983	1,04	0,66
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	12546040	0,73	24749871	0,99	0,73
F14	CR18/B/D7	ICU 18/B/D7	18629402	1,08	21148778	0,85	1,27
F15	CR18/B/D180	ICU 18/B/D180	15140946	0,88	21108311	0,85	1,04
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	7285066	0,42	27652779	1,11	0,38
F20	CR15/B/D7	ICU 12/B/D7	8542752	0,49	28943200	1,16	0,43
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			18843256	1,09	24391825	0,98	1,12

Mitofusin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	16750874,5	cal. mean	40515874		
calibrator			17202799	1,03	43007588	1,06	
E15	CR16/B/D0	ICU 16/B/D0	15823035	0,94	37231612	0,92	1,03
E16	CR16/B/D7	ICU 16/B/D7	15114956	0,90	33408164	0,82	1,09
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			16298950	0,97	38024160	0,94	

Mitofusin 2

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	11020027,5	cal. mean	26252326		
calibrator			11437110	1,04	27774192	1,06	0,98
A1	H1/B	Control 17/B	11963250	1,09	26549130	1,01	1,07
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	23658300	0,90	#####
A4	Sample X/T	Sample X/T	excl.	#####	20212632	0,77	#####
A5	H3/B	Control 3/B	7956990	0,72	22203368	0,85	0,85
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	8956350	0,81	22318120	0,85	0,96
A8	CR4/B/D7	ICU 4/B/D7	5867100	0,53	21786380	0,83	0,64
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	7966530	0,72	24650043	0,94	0,77
A12	CR8/B/D7	ICU 11/B/D7	8601975	0,78	26404414	1,01	0,78
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	6799140	0,62	21712328	0,83	0,75
A16	CR10/B/D7	ICU 13/B/D7	7603155	0,69	23581122	0,90	0,77
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	6459210	0,59	23880825	0,91	0,64
A20	CR11/B/D7	ICU 9/B/D7	6695145	0,61	21025925	0,80	0,76
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			10602945	0,96	24730460	0,94	1,02

Mitofusin 2

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

13475857,5

cal. mean

28899095,5

	calibrator		13024575	0,97	27854350	0,96	1,00
B1	H4/B	Control 12/B	12407265	0,92	28818080	1,00	0,92
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	10045125	0,75	29889932	1,03	0,72
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	10484460	0,78	23959672	0,83	0,94
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	29444064	1,02	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	27218925	0,94	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25503469	0,88	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22796640	0,79	#####
B11	CR9/B/D0	ICU 3/B/D0	10143765	0,75	26684496	0,92	0,82
B12	CR9/B/D7	ICU 3/B/D7	10041615	0,75	27844674	0,96	0,77
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	8912655	0,66	29392896	1,02	0,65
B16	CR12/B/D7	ICU 20/B/D7	11824650	0,88	29492584	1,02	0,86
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	11858940	0,88	28194210	0,98	0,90
B20	CR13/B/D7	ICU 22/B/D7	13396950	0,99	25151087	0,87	1,14
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		13927140	1,03	29943841	1,04	1,00

Mitofusin 2

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C

cal. mean

16597362,5

cal. mean

22029750

	calibrator		15856359	0,96	20662800	0,94	1,02
C1	H15/B	Control 2/B	12694942	0,76	20219460	0,92	0,83
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	14974200	0,90	15286365	0,69	1,30
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	15073822	0,91	14141864	0,64	1,41
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	7036942	0,42	14270928	0,65	0,65
C8	CR3/B/D7	ICU 5/B/D7	9544494	0,58	13239864	0,60	0,96
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	11263974	0,68	18034032	0,82	0,83
C12	CR5/B/D7	ICU 8/B/D7	10391906	0,63	15358512	0,70	0,90
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	12990814	0,78	17619684	0,80	0,98
C16	CR7/B/D7	ICU 6/B/D7	11267608	0,68	17926524	0,81	0,83
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	11900465	0,72	14741100	0,67	1,07
C20	CR20/B/D7	ICU 7/B/D7	13601418	0,82	15625234	0,71	1,16
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		17338366	1,04	23396700	1,06	0,98

Mitofusin 2

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	17590980	cal. mean	27123350		
calibrator			16498328	0,94	21501050	0,79	1,18
D1	H2/B	Control 4/B	13668512	0,78	15463789	0,57	1,36
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	15977016	0,91	15370470	0,57	1,60
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	16747500	0,95	13315368	0,49	1,94
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	16599484	0,94	15176872	0,56	1,69
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	12953996	0,74	12146022	0,45	1,64
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	16019388	0,91	13352832	0,49	1,85
D12	CR2/B/D7	ICU 14/B/D7	15501552	0,88	14137728	0,52	1,69
D13	CR2/B/D180	ICU 14/B/D180	13943776	0,79	11800978	0,44	1,82
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	16267196	0,92	13229508	0,49	1,90
D18	CR17/B/D7	ICU 23/B/D7	12926364	0,73	11261786	0,42	1,77
D19	CR17/B/D180	ICU 23/B/D180	17630096	1,00	10989998	0,41	2,47
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			18683632	1,06	32745650	1,21	0,88

Mitofusin 2

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	17913462,5	cal. mean	13612646,5		
calibrator			18631270	1,04	13056250	0,96	1,08
E1	H12/B	Control 7/B	17375400	0,97	15730101	1,16	0,84
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	11407135	0,64	17836110	1,31	0,49
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	13643160	0,76	17137798	1,26	0,60
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	15734331	0,88	16121115	1,18	0,74
E8	CR6/B/D7	ICU 21/B/D7	11669677	0,65	12428020	0,91	0,71
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	15051280	0,84	12644640	0,93	0,90
E12	CR14/B/D7	ICU 10/B/D7	14287248	0,80	13423200	0,99	0,81
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####
E19	CR19/B/D0	ICU 15/B/D0	15642023	0,87	18376065	1,35	0,65
E20	CR19/B/D7	ICU 15/B/D7	15919605	0,89	18750060	1,38	0,65
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			17195655	0,96	14169043	1,04	0,92

Mitofusin 2

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	18120435	cal. mean	24950732		
calibrator			16884392	0,93	25509639	1,02	0,91
F1	H9/B	Control 9/B	16249224	0,90	32156010	1,29	0,70
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	19760136	1,09	31970444	1,28	0,85
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	18756638	1,04	25759862	1,03	1,00
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	12389640	0,68	24703341	0,99	0,69
F8	CR1/B/D7	ICU 1/B/D7	18136328	1,00	22559436	0,90	1,11
F9	CR1/B/D180	ICU 1/B/D180	12857598	0,71	26013983	1,04	0,68
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	14615074	0,81	24749871	0,99	0,81
F14	CR18/B/D7	ICU 18/B/D7	16344076	0,90	21148778	0,85	1,06
F15	CR18/B/D180	ICU 18/B/D180	18237068	1,01	21108311	0,85	1,19
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	17996810	0,99	27652779	1,11	0,90
F20	CR15/B/D7	ICU 12/B/D7	12544606	0,69	28943200	1,16	0,60
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			19356478	1,07	24391825	0,98	1,09

Mitofusin 2

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	11266359	cal. mean	40515874		
calibrator			10278270	0,91	43007588	1,06	
E15	CR16/B/D0	ICU 16/B/D0	14640816	1,30	37231612	0,92	1,41
E16	CR16/B/D7	ICU 16/B/D7	12435401	1,10	33408164	0,82	1,34
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			12254448	1,09	38024160	0,94	

Tom20

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	4425041,5	cal. mean	25207886		
calibrator			4451682	1,01	26237938	1,04	0,97
A1	H1/B	Control 17/B	17192040	3,89	25525559	1,01	3,84
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	20146728	0,80	#####
A4	Sample X/T	Sample X/T	excl.	#####	19557452	0,78	#####
A5	H3/B	Control 3/B	6596168	1,49	20382535	0,81	1,84
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	11943288	2,70	20062467	0,80	3,39
A8	CR4/B/D7	ICU 4/B/D7	7102452	1,61	19913007	0,79	2,03
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	4212281	0,95	16150422	0,64	1,49
A12	CR8/B/D7	ICU 11/B/D7	8162020	1,84	17634541	0,70	2,64
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	9570422	2,16	17586742	0,70	3,10
A16	CR10/B/D7	ICU 13/B/D7	9561210	2,16	18543098	0,74	2,94
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	8218935	1,86	20294868	0,81	2,31
A20	CR11/B/D7	ICU 9/B/D7	7490956	1,69	21152021	0,84	2,02
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			4398401	0,99	24177834	0,96	1,04

Tom20

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

4671836

cal. mean

28899095,5

	calibrator		4140040	0,89	27854350	0,96	0,92
B1	H4/B	Control 12/B	11238210	2,41	28818080	1,00	2,41
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	7582886	1,62	29889932	1,03	1,57
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	5051560	1,08	23959672	0,83	1,30
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	29444064	1,02	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	27218925	0,94	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25503469	0,88	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22796640	0,79	#####
B11	CR9/B/D0	ICU 3/B/D0	6755616	1,45	26684496	0,92	1,57
B12	CR9/B/D7	ICU 3/B/D7	8607016	1,84	27844674	0,96	1,91
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	7880922	1,69	29392896	1,02	1,66
B16	CR12/B/D7	ICU 20/B/D7	11616288	2,49	29492584	1,02	2,44
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	10723770	2,30	28194210	0,98	2,35
B20	CR13/B/D7	ICU 22/B/D7	11730645	2,51	25151087	0,87	2,89
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		5203632	1,11	29943841	1,04	1,07

Tom20

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 3/ set C		cal. mean	11802124	cal. mean	34079936	
calibrator		12292398	1,04	34138962	1,00	1,04
C1 H15/B	Control 2/B	15812774	1,34	29730870	0,87	1,54
C2 H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3 H16/B	Control 14/B	19865736	1,68	30098518	0,88	1,91
C4 H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5 H17/B	Control 10/B	20022517	1,70	20994664	0,62	2,75
C6 H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7 CR3/B/D0	ICU 5/B/D0	8461645	0,72	22660589	0,66	1,08
C8 CR3/B/D7	ICU 5/B/D7	9581264	0,81	18025875	0,53	1,53
C9 CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10 CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11 CR5/B/D0	ICU 8/B/D0	6937680	0,59	19994035	0,59	1,00
C12 CR5/B/D7	ICU 8/B/D7	8281176	0,70	17391984	0,51	1,37
C13 CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14 CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15 CR7/B/D0	ICU 6/B/D0	11165085	0,95	22085112	0,65	1,46
C16 CR7/B/D7	ICU 6/B/D7	9844808	0,83	18388750	0,54	1,55
C17 CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18 CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19 CR20/B/D0	ICU 7/B/D0	9816420	0,83	24852613	0,73	1,14
C20 CR20/B/D7	ICU 7/B/D7	13241936	1,12	25778115	0,76	1,48
C21 CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22 CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator		11311850	0,96	34020910	1,00	0,96

Tom20

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	8388701		cal. mean	35190284	
calibrator			6974050	0,83	37167787	1,06	0,79
D1	H2/B	Control 4/B	10812209	1,29	30949908	0,88	1,47
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	12141984	1,45	31992571	0,91	1,59
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	17289900	2,06	24406066	0,69	2,97
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	15443580	1,84	22372376	0,64	2,90
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	12351350	1,47	18229138	0,52	2,84
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	13085616	1,56	17679473	0,50	3,10
D12	CR2/B/D7	ICU 14/B/D7	13055130	1,56	18461459	0,52	2,97
D13	CR2/B/D180	ICU 14/B/D180	13513769	1,61	19948365	0,57	2,84
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	11457082	1,37	22814050	0,65	2,11
D18	CR17/B/D7	ICU 23/B/D7	10009395	1,19	24909295	0,71	1,69
D19	CR17/B/D180	ICU 23/B/D180	8937150	1,07	26712850	0,76	1,40
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			9803352	1,17	33212781	0,94	1,24

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	12505153	cal. mean	13612646,5		
calibrator			13380434	1,07	13056250	0,96	1,12
E1	H12/B	Control 7/B	20682106	1,65	15730101	1,16	1,43
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	6727140	0,54	17836110	1,31	0,41
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	12127824	0,97	17137798	1,26	0,77
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	13466267	1,08	16121115	1,18	0,91
E8	CR6/B/D7	ICU 21/B/D7	7086624	0,57	12428020	0,91	0,62
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	14370768	1,15	12644640	0,93	1,24
E12	CR14/B/D7	ICU 10/B/D7	9519427	0,76	13423200	0,99	0,77
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####
E19	CR19/B/D0	ICU 15/B/D0	23031312	1,84	18376065	1,35	1,36
E20	CR19/B/D7	ICU 15/B/D7	18353622	1,47	18750060	1,38	1,07
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			11629872	0,93	14169043	1,04	0,89

Tom20

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	13395526	cal. mean	24950732		
	calibrator		8776741	0,66	25509639	1,02	0,64
F1	H9/B	Control 9/B	23084964	1,72	32156010	1,29	1,34
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	22475025	1,68	31970444	1,28	1,31
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	22393710	1,67	25759862	1,03	1,62
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	12288857	0,92	24703341	0,99	0,93
F8	CR1/B/D7	ICU 1/B/D7	16994474	1,27	22559436	0,90	1,40
F9	CR1/B/D180	ICU 1/B/D180	13326004	0,99	26013983	1,04	0,95
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	7140735	0,53	24749871	0,99	0,54
F14	CR18/B/D7	ICU 18/B/D7	19496475	1,46	21148778	0,85	1,72
F15	CR18/B/D180	ICU 18/B/D180	10199203	0,76	21108311	0,85	0,90
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	13406890	1,00	27652779	1,11	0,90
F20	CR15/B/D7	ICU 12/B/D7	12205260	0,91	28943200	1,16	0,79
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		18014311	1,34	24391825	0,98	1,38

Tom20

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	12471333	cal. mean	40515874		
calibrator			14151445	1,13	43007588	1,06	
E15	CR16/B/D0	ICU 16/B/D0	22650750	1,82	37231612	0,92	1,98
E16	CR16/B/D7	ICU 16/B/D7	21523344	1,73	33408164	0,82	2,09
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			10791221	0,87	38024160	0,94	

Tom40

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	6086490	cal. mean	26252326		
	calibrator		8780756	1,44	27774192	1,06	1,36
A1	H1/B	Control 17/B	13770652	2,26	26549130	1,01	2,24
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	23658300	0,90	#####
A4	Sample X/T	Sample X/T	excl.	#####	20212632	0,77	#####
A5	H3/B	Control 3/B	8353370	1,37	22203368	0,85	1,62
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	9404010	1,55	22318120	0,85	1,82
A8	CR4/B/D7	ICU 4/B/D7	7821932	1,29	21786380	0,83	1,55
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	8265970	1,36	24650043	0,94	1,45
A12	CR8/B/D7	ICU 11/B/D7	8713182	1,43	26404414	1,01	1,42
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	5147124	0,85	21712328	0,83	1,02
A16	CR10/B/D7	ICU 13/B/D7	4873470	0,80	23581122	0,90	0,89
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	2436804	0,40	23880825	0,91	0,44
A20	CR11/B/D7	ICU 9/B/D7	2069080	0,34	21025925	0,80	0,42
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		3392224	0,56	24730460	0,94	0,59

Tom40

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 2/ set B		cal. mean	5630952,5		cal. mean	28899095,5	
calibrator			8553812	1,52	27854350	0,96	1,58
B1	H4/B	Control 12/B	11642840	2,07	28818080	1,00	2,07
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	6467576	1,15	29889932	1,03	1,11
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	8627696	1,53	23959672	0,83	1,85
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	29444064	1,02	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	27218925	0,94	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25503469	0,88	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22796640	0,79	#####
B11	CR9/B/D0	ICU 3/B/D0	7114545	1,26	26684496	0,92	1,37
B12	CR9/B/D7	ICU 3/B/D7	6541366	1,16	27844674	0,96	1,21
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	6018452	1,07	29392896	1,02	1,05
B16	CR12/B/D7	ICU 20/B/D7	7390797	1,31	29492584	1,02	1,29
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	5199656	0,92	28194210	0,98	0,95
B20	CR13/B/D7	ICU 22/B/D7	4362737	0,77	25151087	0,87	0,89
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			2708093	0,48	29943841	1,04	0,46

Tom40

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C

cal. mean

2093057

cal. mean

33101325,5

	calibrator		2885850	1,38	39634150	1,20	1,15
C1	H15/B	Control 2/B	4144545	1,98	41027124	1,24	1,60
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	9156690	4,37	36119890	1,09	4,01
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	13209660	6,31	32155610	0,97	6,50
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	3349980	1,60	31662000	0,96	1,67
C8	CR3/B/D7	ICU 5/B/D7	4827375	2,31	32580420	0,98	2,34
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	4450635	2,13	35399622	1,07	1,99
C12	CR5/B/D7	ICU 8/B/D7	4859595	2,32	28886670	0,87	2,66
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	7071525	3,38	25648220	0,77	4,36
C16	CR7/B/D7	ICU 6/B/D7	9559890	4,57	24801885	0,75	6,10
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	6109246	2,92	23268402	0,70	4,15
C20	CR20/B/D7	ICU 7/B/D7	4242288	2,03	26161994	0,79	2,56
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		1300264	0,62	26568501	0,80	0,77

Tom40

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	2148817,5		cal. mean	37613053,5	
calibrator			2035980	0,95	35832900	0,95	0,99
D1	H2/B	Control 4/B	2928996	1,36	33954570	0,90	1,51
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	7108470	3,31	37204616	0,99	3,34
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	9924885	4,62	27370416	0,73	6,35
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	13499280	6,28	33493056	0,89	7,05
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	10860120	5,05	29291952	0,78	6,49
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	12405105	5,77	26753646	0,71	8,12
D12	CR2/B/D7	ICU 14/B/D7	11848950	5,51	28843383	0,77	7,19
D13	CR2/B/D180	ICU 14/B/D180	16130970	7,51	26981280	0,72	10,46
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	8504100	3,96	28709933	0,76	5,18
D18	CR17/B/D7	ICU 23/B/D7	7105770	3,31	28537335	0,76	4,36
D19	CR17/B/D180	ICU 23/B/D180	5527215	2,57	29401155	0,78	3,29
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			2261655	1,05	39393207	1,05	1,00

Tom40

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	7591203	cal. mean	33459356		
calibrator			6328791	0,83	32367400	0,97	0,86
E1	H12/B	Control 7/B	2810480	0,37	33789330	1,01	0,37
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	405720	0,05	32255844	0,96	0,06
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	1587105	0,21	27424188	0,82	0,26
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	2641995	0,35	29010256	0,87	0,40
E8	CR6/B/D7	ICU 21/B/D7	2254185	0,30	22604584	0,68	0,44
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	6847200	0,90	24246002	0,72	1,24
E12	CR14/B/D7	ICU 10/B/D7	7173450	0,94	26351913	0,79	1,20
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	29617830	0,89	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	26893768	0,80	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	29597504	0,88	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	32157358	0,96	#####
E19	CR19/B/D0	ICU 15/B/D0	13346055	1,76	32525266	0,97	1,81
E20	CR19/B/D7	ICU 15/B/D7	11550240	1,52	38332904	1,15	1,33
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			8853615	1,17	34551312	1,03	1,13

Tom40

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	2152623	cal. mean	37954347		
calibrator			2223522	1,03	26550450	0,70	1,48
F1	H9/B	Control 9/B	2249889	1,05	30236175	0,80	1,31
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	4566156	2,12	28345900	0,75	2,84
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	3759345	1,75	27331912	0,72	2,43
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	3697920	1,72	26053802	0,69	2,50
F8	CR1/B/D7	ICU 1/B/D7	7857900	3,65	27329985	0,72	5,07
F9	CR1/B/D180	ICU 1/B/D180	5578020	2,59	29696880	0,78	3,31
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	5830965	2,71	31147879	0,82	3,30
F14	CR18/B/D7	ICU 18/B/D7	9902025	4,60	34315678	0,90	5,09
F15	CR18/B/D180	ICU 18/B/D180	7274115	3,38	36800640	0,97	3,49
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	1745865	0,81	45890106	1,21	0,67
F20	CR15/B/D7	ICU 12/B/D7	2874105	1,34	39184088	1,03	1,29
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			2081724	0,97	49358244	1,30	0,74

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	5265944	cal. mean	40515874		
calibrator			6566256	1,25	43007588	1,06	
E15	CR16/B/D0	ICU 16/B/D0	10411042	1,98	37231612	0,92	2,15
E16	CR16/B/D7	ICU 16/B/D7	6811542	1,29	33408164	0,82	1,57
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			3965632	0,75	38024160	0,94	

P-DRP1 (Ser 637)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	2992300	cal. mean	32397000		
calibrator			2451386	0,82	34486272	1,06	0,77
A1	H1/B	Control 17/B	4011660	1,34	36056208	1,11	1,20
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	27900540	0,86	#####
A4	Sample X/T	Sample X/T	excl.	#####	28724784	0,89	#####
A5	H3/B	Control 3/B	5213318	1,74	28793205	0,89	1,96
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	2710412	0,91	24152530	0,75	1,21
A8	CR4/B/D7	ICU 4/B/D7	4764220	1,59	22376194	0,69	2,31
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	8360684	2,79	23807098	0,73	3,80
A12	CR8/B/D7	ICU 11/B/D7	8423566	2,82	24407424	0,75	3,74
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	1995296	0,67	23284992	0,72	0,93
A16	CR10/B/D7	ICU 13/B/D7	2249308	0,75	23631488	0,73	1,03
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	2218304	0,74	33965550	1,05	0,71
A20	CR11/B/D7	ICU 9/B/D7	2166554	0,72	32601396	1,01	0,72
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			3533214	1,18	30307728	0,94	1,26

P-DRP1 (Ser 637)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

3298747,5

cal. mean

27865678,5

	calibrator		3525975	1,07	30163392	1,08	0,99
B1	H4/B	Control 12/B	5400450	1,64	27105885	0,97	1,68
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	4220055	1,28	28410290	1,02	1,25
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	4553415	1,38	20371736	0,73	1,89
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	22090275	0,79	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	19395114	0,70	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	20887034	0,75	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	18124046	0,65	#####
B11	CR9/B/D0	ICU 3/B/D0	2910510	0,88	21057888	0,76	1,17
B12	CR9/B/D7	ICU 3/B/D7	4751685	1,44	20873640	0,75	1,92
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	5427585	1,65	26873664	0,96	1,71
B16	CR12/B/D7	ICU 20/B/D7	2650590	0,80	23672428	0,85	0,95
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	2379375	0,72	22632840	0,81	0,89
B20	CR13/B/D7	ICU 22/B/D7	3458925	1,05	20968515	0,75	1,39
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		3071520	0,93	25567965	0,92	1,01

P-DRP1 (Ser 637)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C

cal. mean

3258810***

cal. mean

33101325,5

	calibrator		3258810	1,00	39634150	1,20	0,84
C1	H15/B	Control 2/B	6009696	1,84	41027124	1,24	1,49
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	2742980	0,84	36119890	1,09	0,77
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	1578510	0,48	32155610	0,97	0,50
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	2302785	0,71	31662000	0,96	0,74
C8	CR3/B/D7	ICU 5/B/D7	5476365	1,68	32580420	0,98	1,71
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	6957300	2,13	35399622	1,07	2,00
C12	CR5/B/D7	ICU 8/B/D7	16564050	5,08	28886670	0,87	5,82
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	2771778	0,85	25648220	0,77	1,10
C16	CR7/B/D7	ICU 6/B/D7	1862775	0,57	24801885	0,75	0,76
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	1746585	0,54	23268402	0,70	0,76
C20	CR20/B/D7	ICU 7/B/D7	2664585	0,82	26161994	0,79	1,03
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		967920	0,30	26568501	0,80	0,37

P-DRP1 (Ser 637)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 4/ set D

cal. mean

4142352 ***

cal. mean

37613053,5

	calibrator		4142352	1,00	35832900	0,95	1,05
D1	H2/B	Control 4/B	5240970	1,27	33954570	0,90	1,40
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	4746832	1,15	37204616	0,99	1,16
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	5325075	1,29	27370416	0,73	1,77
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	3935925	0,95	33493056	0,89	1,07
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	3130605	0,76	29291952	0,78	0,97
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	4192408	1,01	26753646	0,71	1,42
D12	CR2/B/D7	ICU 14/B/D7	2185290	0,53	28843383	0,77	0,69
D13	CR2/B/D180	ICU 14/B/D180	3134745	0,76	26981280	0,72	1,05
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	2814165	0,68	28709933	0,76	0,89
D18	CR17/B/D7	ICU 23/B/D7	4329810	1,05	28537335	0,76	1,38
D19	CR17/B/D180	ICU 23/B/D180	10378035	2,51	29401155	0,78	3,21
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		2121936	0,51	39393207	1,05	0,49

P-DRP1 (Ser 637)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	4893312		cal. mean	33459356	
calibrator			4093680	0,84	32367400	0,97	0,86
E1	H12/B	Control 7/B	4672052	0,95	33789330	1,01	0,95
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	8830252	1,80	32255844	0,96	1,87
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	12166020	2,49	27424188	0,82	3,03
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	10031310	2,05	29010256	0,87	2,36
E8	CR6/B/D7	ICU 21/B/D7	13332285	2,72	22604584	0,68	4,03
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	13113628	2,68	24246002	0,72	3,70
E12	CR14/B/D7	ICU 10/B/D7	11918880	2,44	26351913	0,79	3,09
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	CR16/B/D0	ICU 16/B/D0	4153874	0,85	29617830	0,89	0,96
E16	CR16/B/D7	ICU 16/B/D7	2510910	0,51	26893768	0,80	0,64
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E19	CR19/B/D0	ICU 15/B/D0	6654870	1,36	32525266	0,97	1,40
E20	CR19/B/D7	ICU 15/B/D7	5248804	1,07	38332904	1,15	0,94
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			5692944	1,16	34551312	1,03	1,13

P-DRP1 (Ser 637)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	2534847,5		cal. mean	37954347	
	calibrator		2912800	1,15	26550450	0,70	1,64
F1	H9/B	Control 9/B	9266985	3,66	30236175	0,80	4,59
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	3815680	1,51	28345900	0,75	2,02
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	6793160	2,68	27331912	0,72	3,72
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	4377296	1,73	26053802	0,69	2,52
F8	CR1/B/D7	ICU 1/B/D7	6182944	2,44	27329985	0,72	3,39
F9	CR1/B/D180	ICU 1/B/D180	8239275	3,25	29696880	0,78	4,15
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	9189111	3,63	31147879	0,82	4,42
F14	CR18/B/D7	ICU 18/B/D7	2712380	1,07	34315678	0,90	1,18
F15	CR18/B/D180	ICU 18/B/D180	15312976	6,04	36800640	0,97	6,23
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	8114738	3,20	45890106	1,21	2,65
F20	CR15/B/D7	ICU 12/B/D7	21111413	8,33	39184088	1,03	8,07
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		2156895	0,85	49358244	1,30	0,65

DRP1 total (DRP-1 Ser 637 membrane)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	20447491		cal. mean	32397000	
	calibrator		19851813	0,97	34486272	1,06	0,91
A1	H1/B	Control 17/B	16554152	0,81	36056208	1,11	0,73
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	27900540	0,86	#####
A4	Sample X/T	Sample X/T	excl.	#####	28724784	0,89	#####
A5	H3/B	Control 3/B	17454345	0,85	28793205	0,89	0,96
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	19906695	0,97	24152530	0,75	1,31
A8	CR4/B/D7	ICU 4/B/D7	12477344	0,61	22376194	0,69	0,88
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	14319855	0,70	23807098	0,73	0,95
A12	CR8/B/D7	ICU 11/B/D7	17135072	0,84	24407424	0,75	1,11
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	14442795	0,71	23284992	0,72	0,98
A16	CR10/B/D7	ICU 13/B/D7	17411714	0,85	23631488	0,73	1,17
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	17970682	0,88	33965550	1,05	0,84
A20	CR11/B/D7	ICU 9/B/D7	15427454	0,75	32601396	1,01	0,75
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		21043169	1,03	30307728	0,94	1,10

DRP1 total (DRP-1 Ser 637 membrane)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

19760197

cal. mean

27865678,5

	calibrator		19240742	0,97	30163392	1,08	0,90
B1	H4/B	Control 12/B	16654302	0,84	27105885	0,97	0,87
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	13381764	0,68	28410290	1,02	0,66
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	16278416	0,82	20371736	0,73	1,13
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	22090275	0,79	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	19395114	0,70	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	20887034	0,75	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	18124046	0,65	#####
B11	CR9/B/D0	ICU 3/B/D0	15626890	0,79	21057888	0,76	1,05
B12	CR9/B/D7	ICU 3/B/D7	16916776	0,86	20873640	0,75	1,14
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	14959154	0,76	26873664	0,96	0,78
B16	CR12/B/D7	ICU 20/B/D7	17087160	0,86	23672428	0,85	1,02
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	15853596	0,80	22632840	0,81	0,99
B20	CR13/B/D7	ICU 22/B/D7	17871920	0,90	20968515	0,75	1,20
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		20279652	1,03	25567965	0,92	1,12

DRP1 total (DRP-1 Ser 637 membrane)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C
cal. mean
27314950,5
cal. mean
33101325,5

	calibrator		28575180	1,05	39634150	1,20	0,87
C1	H15/B	Control 2/B	28470870	1,04	41027124	1,24	0,84
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	25753024	0,94	36119890	1,09	0,86
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	31206060	1,14	32155610	0,97	1,18
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	18451665	0,68	31662000	0,96	0,71
C8	CR3/B/D7	ICU 5/B/D7	18915705	0,69	32580420	0,98	0,70
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	21493485	0,79	35399622	1,07	0,74
C12	CR5/B/D7	ICU 8/B/D7	22927095	0,84	28886670	0,87	0,96
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	20380635	0,75	25648220	0,77	0,96
C16	CR7/B/D7	ICU 6/B/D7	20218320	0,74	24801885	0,75	0,99
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	22083726	0,81	23268402	0,70	1,15
C20	CR20/B/D7	ICU 7/B/D7	20050704	0,73	26161994	0,79	0,93
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		26054721	0,95	26568501	0,80	1,19

DRP1 total (DRP-1 Ser 637 membrane)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 4/ set D
cal. mean
29612937
cal. mean
37613053,5

	calibrator		32250324	1,09	35832900	0,95	1,14
D1	H2/B	Control 4/B	22447264	0,76	33954570	0,90	0,84
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	32323924	1,09	37204616	0,99	1,10
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	36727826	1,24	27370416	0,73	1,70
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	30768112	1,04	33493056	0,89	1,17
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	23954914	0,81	29291952	0,78	1,04
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	30753668	1,04	26753646	0,71	1,46
D12	CR2/B/D7	ICU 14/B/D7	28668626	0,97	28843383	0,77	1,26
D13	CR2/B/D180	ICU 14/B/D180	26063922	0,88	26981280	0,72	1,23
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	27033878	0,91	28709933	0,76	1,20
D18	CR17/B/D7	ICU 23/B/D7	25482896	0,86	28537335	0,76	1,13
D19	CR17/B/D180	ICU 23/B/D180	28621980	0,97	29401155	0,78	1,24
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		26975550	0,91	39393207	1,05	0,87

DRP1 total (DRP-1 Ser 637 membrane)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	30894397		cal. mean	33459356	
calibrator			32990426	1,07	32367400	0,97	1,10
E1	H12/B	Control 7/B	21289858	0,69	33789330	1,01	0,68
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	18432200	0,60	32255844	0,96	0,62
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	32361885	1,05	27424188	0,82	1,28
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	34403400	1,11	29010256	0,87	1,28
E8	CR6/B/D7	ICU 21/B/D7	27952920	0,90	22604584	0,68	1,34
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	34884315	1,13	24246002	0,72	1,56
E12	CR14/B/D7	ICU 10/B/D7	32954085	1,07	26351913	0,79	1,35
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	CR16/B/D0	ICU 16/B/D0	27389385	0,89	29617830	0,89	1,00
E16	CR16/B/D7	ICU 16/B/D7	28536975	0,92	26893768	0,80	1,15
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E19	CR19/B/D0	ICU 15/B/D0	36932220	1,20	32525266	0,97	1,23
E20	CR19/B/D7	ICU 15/B/D7	37293975	1,21	38332904	1,15	1,05
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			28798368	0,93	34551312	1,03	0,90

DRP1 total (DRP-1 Ser 637 membrane)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	33140997,5		cal. mean	37954347	
calibrator			34865043	1,05	26550450	0,70	1,50
F1	H9/B	Control 9/B	32431635	0,98	30236175	0,80	1,23
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	39592293	1,19	28345900	0,75	1,60
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	40974030	1,24	27331912	0,72	1,72
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	24593760	0,74	26053802	0,69	1,08
F8	CR1/B/D7	ICU 1/B/D7	34650810	1,05	27329985	0,72	1,45
F9	CR1/B/D180	ICU 1/B/D180	31007616	0,94	29696880	0,78	1,20
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	27601008	0,83	31147879	0,82	1,01
F14	CR18/B/D7	ICU 18/B/D7	30641585	0,92	34315678	0,90	1,02
F15	CR18/B/D180	ICU 18/B/D180	36203610	1,09	36800640	0,97	1,13
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	34813471	1,05	45890106	1,21	0,87
F20	CR15/B/D7	ICU 12/B/D7	28861470	0,87	39184088	1,03	0,84
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			31416952	0,95	49358244	1,30	0,73

P-DRP1(Ser 637) vs. DRP1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
GEL 1/ set A		cal. mean	2992300	cal. mean	20447491		
calibrator			2451386	0,82	19851813	0,97	0,84
A1	H1/B	Control 17/B	4011660	1,34	16554152	0,81	1,66
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	14963011	0,73	#####
A4	Sample X/T	Sample X/T	excl.	#####	6715595	0,33	#####
A5	H3/B	Control 3/B	5213318	1,74	17454345	0,85	2,04
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	2710412	0,91	19906695	0,97	0,93
A8	CR4/B/D7	ICU 4/B/D7	4764220	1,59	12477344	0,61	2,61
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	8360684	2,79	14319855	0,70	3,99
A12	CR8/B/D7	ICU 11/B/D7	8423566	2,82	17135072	0,84	3,36
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	1995296	0,67	14442795	0,71	0,94
A16	CR10/B/D7	ICU 13/B/D7	2249308	0,75	17411714	0,85	0,88
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	2218304	0,74	17970682	0,88	0,84
A20	CR11/B/D7	ICU 9/B/D7	2166554	0,72	15427454	0,75	0,96
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			3533214	1,18	21043169	1,03	1,15

P-DRP1(Ser 637) vs. DRP1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
-------	--------------------	-----------------	----------------------------------	--------------	-------------------------------	-------------

GEL 2/ set B

cal. mean

3298747,5

cal. mean

19760197

	calibrator		3525975	1,07	19240742	0,97	1,10
B1	H4/B	Control 12/B	5400450	1,64	16654302	0,84	1,94
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	4220055	1,28	13381764	0,68	1,89
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	4553415	1,38	16278416	0,82	1,68
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	14014774	0,71	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	10927346	0,55	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	13455092	0,68	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	10346826	0,52	#####
B11	CR9/B/D0	ICU 3/B/D0	2910510	0,88	15626890	0,79	1,12
B12	CR9/B/D7	ICU 3/B/D7	4751685	1,44	16916776	0,86	1,68
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	5427585	1,65	14959154	0,76	2,17
B16	CR12/B/D7	ICU 20/B/D7	2650590	0,80	17087160	0,86	0,93
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	2379375	0,72	15853596	0,80	0,90
B20	CR13/B/D7	ICU 22/B/D7	3458925	1,05	17871920	0,90	1,16
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		3071520	0,93	20279652	1,03	0,91

P-DRP1(Ser 637) vs. DRP1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
-------	--------------------	-----------------	----------------------------------	--------------	-------------------------------	-------------

GEL 3/ set C

cal. mean 3258810 ***

cal. mean 27314950,5

	calibrator		3258810	1,00	28575180	1,05	0,96
C1	H15/B	Control 2/B	6009696	1,84	28470870	1,04	1,77
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	2742980	0,84	25753024	0,94	0,89
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	1578510	0,48	31206060	1,14	0,42
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	2302785	0,71	18451665	0,68	1,05
C8	CR3/B/D7	ICU 5/B/D7	5476365	1,68	18915705	0,69	2,43
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	6957300	2,13	21493485	0,79	2,71
C12	CR5/B/D7	ICU 8/B/D7	16564050	5,08	22927095	0,84	6,06
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	2771778	0,85	20380635	0,75	1,14
C16	CR7/B/D7	ICU 6/B/D7	1862775	0,57	20218320	0,74	0,77
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	1746585	0,54	22083726	0,81	0,66
C20	CR20/B/D7	ICU 7/B/D7	2664585	0,82	20050704	0,73	1,11
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		967920	0,30	26054721	0,95	0,31

P-DRP1(Ser 637) vs. DRP1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
-------	--------------------	-----------------	----------------------------------	--------------	-------------------------------	-------------

GEL 4/ set D

cal. mean

4142352 ***

cal. mean

29612937

	calibrator		4142352	1,00	32250324	1,09	0,92
D1	H2/B	Control 4/B	5240970	1,27	22447264	0,76	1,67
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	4746832	1,15	32323924	1,09	1,05
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	5325075	1,29	36727826	1,24	1,04
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	3935925	0,95	30768112	1,04	0,91
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	3130605	0,76	23954914	0,81	0,93
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	4192408	1,01	30753668	1,04	0,97
D12	CR2/B/D7	ICU 14/B/D7	2185290	0,53	28668626	0,97	0,54
D13	CR2/B/D180	ICU 14/B/D180	3134745	0,76	26063922	0,88	0,86
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	2814165	0,68	27033878	0,91	0,74
D18	CR17/B/D7	ICU 23/B/D7	4329810	1,05	25482896	0,86	1,21
D19	CR17/B/D180	ICU 23/B/D180	10378035	2,51	28621980	0,97	2,59
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		2121936	0,51	26975550	0,91	0,56

P-DRP1(Ser 637) vs. DRP1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
GEL 5/ set E		cal. mean	4893312		cal. mean	30894397	
calibrator			4093680	0,84	32990426	1,07	0,78
E1	H12/B	Control 7/B	4672052	0,95	21289858	0,69	1,39
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	8830252	1,80	18432200	0,60	3,02
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	12166020	2,49	32361885	1,05	2,37
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	10031310	2,05	34403400	1,11	1,84
E8	CR6/B/D7	ICU 21/B/D7	13332285	2,72	27952920	0,90	3,01
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	13113628	2,68	34884315	1,13	2,37
E12	CR14/B/D7	ICU 10/B/D7	11918880	2,44	32954085	1,07	2,28
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	CR16/B/D0	ICU 16/B/D0	4153874	0,85	27389385	0,89	0,96
E16	CR16/B/D7	ICU 16/B/D7	2510910	0,51	28536975	0,92	0,56
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E19	CR19/B/D0	ICU 15/B/D0	6654870	1,36	36932220	1,20	1,14
E20	CR19/B/D7	ICU 15/B/D7	5248804	1,07	37293975	1,21	0,89
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			5692944	1,16	28798368	0,93	1,25

P-DRP1(Ser 637) vs. DRP1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
-------	--------------------	-----------------	----------------------------------	--------------	-------------------------------	-------------

GEL 6/ set F

cal. mean

2534847,5

cal. mean

33140997,5

	calibrator		2912800	1,15	34865043	1,05	1,09
F1	H9/B	Control 9/B	9266985	3,66	32431635	0,98	3,74
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	3815680	1,51	39592293	1,19	1,26
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	6793160	2,68	40974030	1,24	2,17
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	4377296	1,73	24593760	0,74	2,33
F8	CR1/B/D7	ICU 1/B/D7	6182944	2,44	34650810	1,05	2,33
F9	CR1/B/D180	ICU 1/B/D180	8239275	3,25	31007616	0,94	3,47
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	9189111	3,63	27601008	0,83	4,35
F14	CR18/B/D7	ICU 18/B/D7	2712380	1,07	30641585	0,92	1,16
F15	CR18/B/D180	ICU 18/B/D180	15312976	6,04	36203610	1,09	5,53
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	8114738	3,20	34813471	1,05	3,05
F20	CR15/B/D7	ICU 12/B/D7	21111413	8,33	28861470	0,87	9,56
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		2156895	0,85	31416952	0,95	0,90

P-DRP1 (Ser616)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	13155333	cal. mean	10427420		
	calibrator		13478092	1,02	11665917	1,12	0,92
A1	H1/B	Control 17/B	1860378	0,14	14751528	1,41	0,10
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	11519366	1,10	#####
A4	Sample X/T	Sample X/T	excl.	#####	9972366	0,96	#####
A5	H3/B	Control 3/B	5922638	0,45	12429150	1,19	0,38
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	1236102	0,09	13267540	1,27	0,07
A8	CR4/B/D7	ICU 4/B/D7	1361094	0,10	14906802	1,43	0,07
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	5276788	0,40	12271770	1,18	0,34
A12	CR8/B/D7	ICU 11/B/D7	5430960	0,41	13834497	1,33	0,31
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	2116845	0,16	16572059	1,59	0,10
A16	CR10/B/D7	ICU 13/B/D7	3168836	0,24	17226359	1,65	0,15
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	1480469	0,11	16502540	1,58	0,07
A20	CR11/B/D7	ICU 9/B/D7	1506204	0,11	16566748	1,59	0,07
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		12832574	0,98	9188923	0,88	1,11

P-DRP1 (Ser616)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

15021821

cal. mean

13620209

	calibrator		15325498	1,02	15103226	1,11	0,92
B1	H4/B	Control 12/B	2279129	0,15	24508110	1,80	0,08
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	1248660	0,08	20180108	1,48	0,06
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	1361272	0,09	15163650	1,11	0,08
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	15430378	1,13	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	13459094	0,99	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	14559168	1,07	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	12338764	0,91	#####
B11	CR9/B/D0	ICU 3/B/D0	3198150	0,21	15477344	1,14	0,19
B12	CR9/B/D7	ICU 3/B/D7	4367838	0,29	19214250	1,41	0,21
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	2361180	0,16	19783015	1,45	0,11
B16	CR12/B/D7	ICU 20/B/D7	3132554	0,21	20225326	1,48	0,14
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	1530225	0,10	16459904	1,21	0,08
B20	CR13/B/D7	ICU 22/B/D7	1226958	0,08	15185244	1,11	0,07
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		14718144	0,98	12137192	0,89	1,10

P-DRP1 (Ser616)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C

cal. mean

24661704 ***

cal. mean

29009047,5

	calibrator		24661704	1,00	31268520	1,08	0,93
C1	H15/B	Control 2/B	3797892	0,15	32658660	1,13	0,14
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	1220288	0,05	30500228	1,05	0,05
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	1466572	0,06	26647155	0,92	0,06
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	471086	0,02	23008587	0,79	0,02
C8	CR3/B/D7	ICU 5/B/D7	887386	0,04	20825559	0,72	0,05
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	1312058	0,05	24818544	0,86	0,06
C12	CR5/B/D7	ICU 8/B/D7	1467354	0,06	21221460	0,73	0,08
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	804126	0,03	31349850	1,08	0,03
C16	CR7/B/D7	ICU 6/B/D7	779562	0,03	30200625	1,04	0,03
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	2678718	0,11	29807955	1,03	0,11
C20	CR20/B/D7	ICU 7/B/D7	5670060	0,23	32177728	1,11	0,21
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		12151020	0,49	26749575	0,92	0,53

P-DRP1 (Ser616)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	24102166		cal. mean	29637461,5	
calibrator			25517316	1,06	32687753	1,10	0,96
D1	H2/B	Control 4/B	3355880	0,14	25631280	0,86	0,16
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	2171444	0,09	31143112	1,05	0,09
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	6775340	0,28	25175295	0,85	0,33
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	2748064	0,11	25479814	0,86	0,13
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	1118304	0,05	20932520	0,71	0,07
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	2075172	0,09	30241344	1,02	0,08
D12	CR2/B/D7	ICU 14/B/D7	1012176	0,04	25133085	0,85	0,05
D13	CR2/B/D180	ICU 14/B/D180	806300	0,03	24795792	0,84	0,04
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	942656	0,04	26724960	0,90	0,04
D18	CR17/B/D7	ICU 23/B/D7	651024	0,03	26946765	0,91	0,03
D19	CR17/B/D180	ICU 23/B/D180	2268332	0,09	30084089	1,02	0,09
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			22687016	0,94	26587170	0,90	1,05

P-DRP1 (Ser616)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 5/ set E

cal. mean

14855670

cal. mean

10923817,5

	calibrator		15512445	1,04	12237750	1,12	0,93
E1	H12/B	Control 7/B	1026045	0,07	12498750	1,14	0,06
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	478890	0,03	19614149	1,80	0,02
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	3126420	0,21	17353592	1,59	0,13
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	4776120	0,32	16839405	1,54	0,21
E8	CR6/B/D7	ICU 21/B/D7	2857770	0,19	16236244	1,49	0,13
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	3490875	0,23	19310702	1,77	0,13
E12	CR14/B/D7	ICU 10/B/D7	3847005	0,26	19300910	1,77	0,15
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	CR16/B/D0	ICU 16/B/D0	3723255	0,25	14058180	1,29	0,19
E16	CR16/B/D7	ICU 16/B/D7	3113955	0,21	13015838	1,19	0,18
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E19	CR19/B/D0	ICU 15/B/D0	4733730	0,32	15185278	1,39	0,23
E20	CR19/B/D7	ICU 15/B/D7	2046780	0,14	14273055	1,31	0,11
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		14198895	0,96	9609885	0,88	1,09

P-DRP1 (Ser616)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 6/ set F

cal. mean

19134292,5

cal. mean

14481045

	calibrator		20892240	1,09	11853990	0,82	1,33
F1	H9/B	Control 9/B	2622285	0,14	19641015	1,36	0,10
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	4891275	0,26	20068098	1,39	0,18
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	4027770	0,21	17405802	1,20	0,18
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	2574405	0,13	15359580	1,06	0,13
F8	CR1/B/D7	ICU 1/B/D7	5365350	0,28	16966201	1,17	0,24
F9	CR1/B/D180	ICU 1/B/D180	3132225	0,16	16526986	1,14	0,14
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	3284910	0,17	16144528	1,11	0,15
F14	CR18/B/D7	ICU 18/B/D7	5788350	0,30	18203006	1,26	0,24
F15	CR18/B/D180	ICU 18/B/D180	5519340	0,29	20387890	1,41	0,20
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	3096720	0,16	18896984	1,30	0,12
F20	CR15/B/D7	ICU 12/B/D7	2158200	0,11	18737415	1,29	0,09
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		17376345	0,91	17108100	1,18	0,77

DRP1 total (P-DRP1 Ser616 membrane)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 1/ set A
cal. mean
20490527
cal. mean
10427420

	calibrator		21299080	1,04	11665917	1,12	0,93
A1	H1/B	Control 17/B	20846716	1,02	14751528	1,41	0,72
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	11519366	1,10	#####
A4	Sample X/T	Sample X/T	excl.	#####	9972366	0,96	#####
A5	H3/B	Control 3/B	18109476	0,88	12429150	1,19	0,74
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	17183540	0,84	13267540	1,27	0,66
A8	CR4/B/D7	ICU 4/B/D7	15410384	0,75	14906802	1,43	0,53
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	19115756	0,93	12271770	1,18	0,79
A12	CR8/B/D7	ICU 11/B/D7	17281968	0,84	13834497	1,33	0,64
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	14008324	0,68	16572059	1,59	0,43
A16	CR10/B/D7	ICU 13/B/D7	18447044	0,90	17226359	1,65	0,54
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	18058348	0,88	16502540	1,58	0,56
A20	CR11/B/D7	ICU 9/B/D7	15125308	0,74	16566748	1,59	0,46
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		19681974	0,96	9188923	0,88	1,09

DRP1 total (P-DRP1 Ser616 membrane)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 2/ set B		cal. mean	22053389,5		cal. mean	13620209	
calibrator			20342028	0,92	15103226	1,11	0,83
B1	H4/B	Control 12/B	21041320	0,95	24508110	1,80	0,53
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	19388310	0,88	20180108	1,48	0,59
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	18330556	0,83	15163650	1,11	0,75
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	15430378	1,13	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	13459094	0,99	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	14559168	1,07	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	12338764	0,91	#####
B11	CR9/B/D0	ICU 3/B/D0	20546084	0,93	15477344	1,14	0,82
B12	CR9/B/D7	ICU 3/B/D7	21659882	0,98	19214250	1,41	0,70
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	20499532	0,93	19783015	1,45	0,64
B16	CR12/B/D7	ICU 20/B/D7	22765354	1,03	20225326	1,48	0,70
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	20832075	0,94	16459904	1,21	0,78
B20	CR13/B/D7	ICU 22/B/D7	22226372	1,01	15185244	1,11	0,90
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			23764751	1,08	12137192	0,89	1,21

DRP1 total (P-DRP1 Ser616 membrane)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 3/ set C		cal. mean	33811909		cal. mean	29009047,5	
calibrator			33888062	1,00	31268520	1,08	0,93
C1	H15/B	Control 2/B	36855246	1,09	32658660	1,13	0,97
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	32162452	0,95	30500228	1,05	0,90
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	37903540	1,12	26647155	0,92	1,22
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	19575254	0,58	23008587	0,79	0,73
C8	CR3/B/D7	ICU 5/B/D7	21974338	0,65	20825559	0,72	0,91
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	33517302	0,99	24818544	0,86	1,16
C12	CR5/B/D7	ICU 8/B/D7	41777016	1,24	21221460	0,73	1,69
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	36343006	1,07	31349850	1,08	0,99
C16	CR7/B/D7	ICU 6/B/D7	30973320	0,92	30200625	1,04	0,88
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	41217886	1,22	29807955	1,03	1,19
C20	CR20/B/D7	ICU 7/B/D7	41811378	1,24	32177728	1,11	1,11
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			33735756	1,00	26749575	0,92	1,08

DRP1 total (P-DRP1 Ser616 membrane)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	29376405		cal. mean	29637461,5	
calibrator			31913742	1,09	32687753	1,10	0,98
D1	H2/B	Control 4/B	32334826	1,10	25631280	0,86	1,27
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	44026324	1,50	31143112	1,05	1,43
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	44171512	1,50	25175295	0,85	1,77
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	39353184	1,34	25479814	0,86	1,56
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	19317984	0,66	20932520	0,71	0,93
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	34569920	1,18	30241344	1,02	1,15
D12	CR2/B/D7	ICU 14/B/D7	30520172	1,04	25133085	0,85	1,23
D13	CR2/B/D180	ICU 14/B/D180	29586878	1,01	24795792	0,84	1,20
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	34750194	1,18	26724960	0,90	1,31
D18	CR17/B/D7	ICU 23/B/D7	30497808	1,04	26946765	0,91	1,14
D19	CR17/B/D180	ICU 23/B/D180	39151106	1,33	30084089	1,02	1,31
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			26839068	0,91	26587170	0,90	1,02

DRP1 total (P-DRP1 Ser616 membrane)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 5/ set E
cal. mean
23540647,5
cal. mean
10923817,5

	calibrator		25810290	1,10	12237750	1,12	0,98
E1	H12/B	Control 7/B	21310380	0,91	12498750	1,14	0,79
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	17777655	0,76	19614149	1,80	0,42
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	32708205	1,39	17353592	1,59	0,87
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	30594600	1,30	16839405	1,54	0,84
E8	CR6/B/D7	ICU 21/B/D7	29838240	1,27	16236244	1,49	0,85
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	36919890	1,57	19310702	1,77	0,89
E12	CR14/B/D7	ICU 10/B/D7	40245075	1,71	19300910	1,77	0,97
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	CR16/B/D0	ICU 16/B/D0	39698092	1,69	14058180	1,29	1,31
E16	CR16/B/D7	ICU 16/B/D7	31497255	1,34	13015838	1,19	1,12
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E19	CR19/B/D0	ICU 15/B/D0	40211595	1,71	15185278	1,39	1,23
E20	CR19/B/D7	ICU 15/B/D7	42000570	1,78	14273055	1,31	1,37
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		21271005	0,90	9609885	0,88	1,03

DRP1 total (P-DRP1 Ser616 membrane)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 6/ set F
cal. mean
29816525,5
cal. mean
14481045

	calibrator		29337650	0,98	11853990	0,82	1,20
F1	H9/B	Control 9/B	30040944	1,01	19641015	1,36	0,74
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	49933506	1,67	20068098	1,39	1,21
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	40714872	1,37	17405802	1,20	1,14
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	23795892	0,80	15359580	1,06	0,75
F8	CR1/B/D7	ICU 1/B/D7	33495544	1,12	16966201	1,17	0,96
F9	CR1/B/D180	ICU 1/B/D180	30412808	1,02	16526986	1,14	0,89
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	17583776	0,59	16144528	1,11	0,53
F14	CR18/B/D7	ICU 18/B/D7	21266996	0,71	18203006	1,26	0,57
F15	CR18/B/D180	ICU 18/B/D180	38725146	1,30	20387890	1,41	0,92
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	26311494	0,88	18896984	1,30	0,68
F20	CR15/B/D7	ICU 12/B/D7	22636504	0,76	18737415	1,29	0,59
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		30295401	1,02	17108100	1,18	0,86

P-DRP1 (S616) vs. DRP1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
GEL 1/ set A		cal. mean	13155333	cal. mean	20490527		
	calibrator		13478092	1,02	21299080	1,04	0,99
A1	H1/B	Control 17/B	1860378	0,14	20846716	1,02	0,14
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	16393960	0,80	#####
A4	Sample X/T	Sample X/T	excl.	#####	5640624	0,28	#####
A5	H3/B	Control 3/B	5922638	0,45	18109476	0,88	0,51
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	1236102	0,09	17183540	0,84	0,11
A8	CR4/B/D7	ICU 4/B/D7	1361094	0,10	15410384	0,75	0,14
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	5276788	0,40	19115756	0,93	0,43
A12	CR8/B/D7	ICU 11/B/D7	5430960	0,41	17281968	0,84	0,49
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	2116845	0,16	14008324	0,68	0,24
A16	CR10/B/D7	ICU 13/B/D7	3168836	0,24	18447044	0,90	0,27
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	1480469	0,11	18058348	0,88	0,13
A20	CR11/B/D7	ICU 9/B/D7	1506204	0,11	15125308	0,74	0,16
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		12832574	0,98	19681974	0,96	1,02

P-DRP1 (S616) vs. DRP1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
-------	--------------------	-----------------	----------------------------------	--------------	-------------------------------	-------------

GEL 2/ set B

cal. mean

15021821

cal. mean

22053389,5

	calibrator		15325498	1,02	20342028	0,92	1,11
B1	H4/B	Control 12/B	2279129	0,15	21041320	0,95	0,16
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	1248660	0,08	19388310	0,88	0,09
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	1361272	0,09	18330556	0,83	0,11
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	18612520	0,84	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	13257016	0,60	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	19560120	0,89	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	15995350	0,73	#####
B11	CR9/B/D0	ICU 3/B/D0	3198150	0,21	20546084	0,93	0,23
B12	CR9/B/D7	ICU 3/B/D7	4367838	0,29	21659882	0,98	0,30
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	2361180	0,16	20499532	0,93	0,17
B16	CR12/B/D7	ICU 20/B/D7	3132554	0,21	22765354	1,03	0,20
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	1530225	0,10	20832075	0,94	0,11
B20	CR13/B/D7	ICU 22/B/D7	1226958	0,08	22226372	1,01	0,08
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		14718144	0,98	23764751	1,08	0,91

P-DRP1 (S616) vs. DRP1

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
-------	--------------------	-----------------	----------------------------------	--------------	-------------------------------	-------------

GEL 3/ set C

cal. mean 24661704 ***

cal. mean 33811909

	calibrator		24661704	1,00	33888062	1,00	1,00
C1	H15/B	Control 2/B	3797892	0,15	36855246	1,09	0,14
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	1220288	0,05	32162452	0,95	0,05
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	1466572	0,06	37903540	1,12	0,05
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	471086	0,02	19575254	0,58	0,03
C8	CR3/B/D7	ICU 5/B/D7	887386	0,04	21974338	0,65	0,06
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	1312058	0,05	33517302	0,99	0,05
C12	CR5/B/D7	ICU 8/B/D7	1467354	0,06	41777016	1,24	0,05
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	804126	0,03	36343006	1,07	0,03
C16	CR7/B/D7	ICU 6/B/D7	779562	0,03	30973320	0,92	0,03
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	2678718	0,11	41217886	1,22	0,09
C20	CR20/B/D7	ICU 7/B/D7	5670060	0,23	41811378	1,24	0,19
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		12151020	0,49	33735756	1,00	0,49

P-DRP1 (S616) vs. DRP1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
GEL 4/ set D		cal. mean	24102166		cal. mean	29637461,5	
calibrator			25517316	1,06	32687753	1,10	0,96
D1	H2/B	Control 4/B	3355880	0,14	25631280	0,86	0,16
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	2171444	0,09	31143112	1,05	0,09
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	6775340	0,28	25175295	0,85	0,33
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	2748064	0,11	25479814	0,86	0,13
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	1118304	0,05	20932520	0,71	0,07
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	2075172	0,09	30241344	1,02	0,08
D12	CR2/B/D7	ICU 14/B/D7	1012176	0,04	25133085	0,85	0,05
D13	CR2/B/D180	ICU 14/B/D180	806300	0,03	24795792	0,84	0,04
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	942656	0,04	26724960	0,90	0,04
D18	CR17/B/D7	ICU 23/B/D7	651024	0,03	26946765	0,91	0,03
D19	CR17/B/D180	ICU 23/B/D180	2268332	0,09	30084089	1,02	0,09
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			22687016	0,94	26587170	0,90	1,05

P-DRP1 (S616) vs. DRP1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
GEL 5/ set E		cal. mean	14855670	cal. mean	23540647,5		
	calibrator		15512445	1,04	25810290	1,10	0,95
E1	H12/B	Control 7/B	1026045	0,07	21310380	0,91	0,08
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	478890	0,03	17777655	0,76	0,04
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	3126420	0,21	32708205	1,39	0,15
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	4776120	0,32	30594600	1,30	0,25
E8	CR6/B/D7	ICU 21/B/D7	2857770	0,19	29838240	1,27	0,15
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	3490875	0,23	36919890	1,57	0,15
E12	CR14/B/D7	ICU 10/B/D7	3847005	0,26	40245075	1,71	0,15
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	CR16/B/D0	ICU 16/B/D0	3723255	0,25	39698092	1,69	0,15
E16	CR16/B/D7	ICU 16/B/D7	3113955	0,21	31497255	1,34	0,16
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E19	CR19/B/D0	ICU 15/B/D0	4733730	0,32	40211595	1,71	0,19
E20	CR19/B/D7	ICU 15/B/D7	2046780	0,14	42000570	1,78	0,08
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		14198895	0,96	21271005	0,90	1,06

P-DRP1 (S616) vs. DRP1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	DRP1 density	DRP1 normalized to calibrator	P-DRP1/DRP1
GEL 6/ set F		cal. mean	19134292,5	cal. mean	29816525,5		
	calibrator		20892240	1,09	29337650	0,98	1,11
F1	H9/B	Control 9/B	2622285	0,14	30040944	1,01	0,14
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	4891275	0,26	49933506	1,67	0,15
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	4027770	0,21	40714872	1,37	0,15
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	2574405	0,13	23795892	0,80	0,17
F8	CR1/B/D7	ICU 1/B/D7	5365350	0,28	33495544	1,12	0,25
F9	CR1/B/D180	ICU 1/B/D180	3132225	0,16	30412808	1,02	0,16
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	3284910	0,17	17583776	0,59	0,29
F14	CR18/B/D7	ICU 18/B/D7	5788350	0,30	21266996	0,71	0,42
F15	CR18/B/D180	ICU 18/B/D180	5519340	0,29	38725146	1,30	0,22
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	3096720	0,16	26311494	0,88	0,18
F20	CR15/B/D7	ICU 12/B/D7	2158200	0,11	22636504	0,76	0,15
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		17376345	0,91	30295401	1,02	0,89

OPA1-L (upper band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 1/ set A
cal. mean
6214612,5
cal. mean
41570172

	calibrator		7497630	1,21	42065400	1,01	1,19
A1	H1/B	Control 17/B	7544340	1,21	42967628	1,03	1,17
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	36839820	0,89	#####
A4	Sample X/T	Sample X/T	excl.	#####	32757912	0,79	#####
A5	H3/B	Control 3/B	7926660	1,28	35733870	0,86	1,48
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	5397570	0,87	28368765	0,68	1,27
A8	CR4/B/D7	ICU 4/B/D7	5267475	0,85	30026880	0,72	1,17
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	5111010	0,82	26405840	0,64	1,29
A12	CR8/B/D7	ICU 11/B/D7	5631615	0,91	28186268	0,68	1,34
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	5443843	0,88	32125544	0,77	1,13
A16	CR10/B/D7	ICU 13/B/D7	5619645	0,90	37040626	0,89	1,01
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	5215408	0,84	34330306	0,83	1,02
A20	CR11/B/D7	ICU 9/B/D7	6976578	1,12	33386188	0,80	1,40
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		4931595	0,79	41074944	0,99	0,80

OPA1-L (upper band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

6667057

cal. mean

39443078,5

	calibrator		6068835	0,91	31672407	0,80	1,13
B1	H4/B	Control 12/B	5838840	0,88	34761556	0,88	0,99
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	5453280	0,82	31297500	0,79	1,03
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	5134680	0,77	26830935	0,68	1,13
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	25290990	0,64	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	23090160	0,59	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25183170	0,64	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22048708	0,56	#####
B11	CR9/B/D0	ICU 3/B/D0	3499245	0,52	28875120	0,73	0,72
B12	CR9/B/D7	ICU 3/B/D7	3088710	0,46	32768164	0,83	0,56
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	3439484	0,52	32436270	0,82	0,63
B16	CR12/B/D7	ICU 20/B/D7	4561744	0,68	33377085	0,85	0,81
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	5225616	0,78	35632337	0,90	0,87
B20	CR13/B/D7	ICU 22/B/D7	5349544	0,80	35268420	0,89	0,90
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		7265279	1,09	47213750	1,20	0,91

OPA1-L (upper band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C
cal. mean
12017075
cal. mean
22029750

	calibrator		11766732	0,98	20662800	0,94	1,04
C1	H15/B	Control 2/B	11093084	0,92	20219460	0,92	1,01
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	11870276	0,99	15286365	0,69	1,42
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	11324460	0,94	14141864	0,64	1,47
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	7537652	0,63	14270928	0,65	0,97
C8	CR3/B/D7	ICU 5/B/D7	9888942	0,82	13239864	0,60	1,37
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	12336970	1,03	18034032	0,82	1,25
C12	CR5/B/D7	ICU 8/B/D7	10709352	0,89	15358512	0,70	1,28
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	12870968	1,07	17619684	0,80	1,34
C16	CR7/B/D7	ICU 6/B/D7	10077036	0,84	17926524	0,81	1,03
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	12286003	1,02	14741100	0,67	1,53
C20	CR20/B/D7	ICU 7/B/D7	14452924	1,20	15625234	0,71	1,70
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		12267418	1,02	23396700	1,06	0,96

OPA1-L (upper band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 4/ set D
cal. mean
9487102,5
cal. mean
27123350

	calibrator		10049400	1,06	21501050	0,79	1,34
D1	H2/B	Control 4/B	12756555	1,34	15463789	0,57	2,36
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	5276070	0,56	15370470	0,57	0,98
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	4978980	0,52	13315368	0,49	1,07
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	5876325	0,62	15176872	0,56	1,11
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	7436700	0,78	12146022	0,45	1,75
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	8344395	0,88	13352832	0,49	1,79
D12	CR2/B/D7	ICU 14/B/D7	7844400	0,83	14137728	0,52	1,59
D13	CR2/B/D180	ICU 14/B/D180	8646750	0,91	11800978	0,44	2,09
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	13831470	1,46	13229508	0,49	2,99
D18	CR17/B/D7	ICU 23/B/D7	11146500	1,17	11261786	0,42	2,83
D19	CR17/B/D180	ICU 23/B/D180	13643550	1,44	10989998	0,41	3,55
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		8924805	0,94	32745650	1,21	0,78

OPA1-L (upper band)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	6622605,5	cal. mean	13612646,5		
calibrator			6786894	1,02	13056250	0,96	1,07
E1	H12/B	Control 7/B	9449632	1,43	15730101	1,16	1,23
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	6779204	1,02	17836110	1,31	0,78
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	7881524	1,19	17137798	1,26	0,95
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	9265439	1,40	16121115	1,18	1,18
E8	CR6/B/D7	ICU 21/B/D7	5076990	0,77	12428020	0,91	0,84
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	6513435	0,98	12644640	0,93	1,06
E12	CR14/B/D7	ICU 10/B/D7	7567282	1,14	13423200	0,99	1,16
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####
E19	CR19/B/D0	ICU 15/B/D0	8358292	1,26	18376065	1,35	0,93
E20	CR19/B/D7	ICU 15/B/D7	7947588	1,20	18750060	1,38	0,87
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			6458317	0,98	14169043	1,04	0,94

OPA1-L (upper band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 6/ set F
cal. mean
9168219,5
cal. mean
24950732

	calibrator		9497525	1,04	25509639	1,02	1,01
F1	H9/B	Control 9/B	12429009	1,36	32156010	1,29	1,05
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	10357650	1,13	31970444	1,28	0,88
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	8423622	0,92	25759862	1,03	0,89
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	7570035	0,83	24703341	0,99	0,83
F8	CR1/B/D7	ICU 1/B/D7	5860538	0,64	22559436	0,90	0,71
F9	CR1/B/D180	ICU 1/B/D180	6842224	0,75	26013983	1,04	0,72
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	9948240	1,09	24749871	0,99	1,09
F14	CR18/B/D7	ICU 18/B/D7	10238950	1,12	21148778	0,85	1,32
F15	CR18/B/D180	ICU 18/B/D180	9353295	1,02	21108311	0,85	1,21
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	14508383	1,58	27652779	1,11	1,43
F20	CR15/B/D7	ICU 12/B/D7	11674457	1,27	28943200	1,16	1,10
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		8838914	0,96	24391825	0,98	0,99

OPA1-L (upper band)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	4300663,5	cal. mean	15623570		
calibrator			4627422	1,08	17691250	1,13	
E15	CR16/B/D0	ICU 16/B/D0	4714126	1,10	14143732	0,91	1,21
E16	CR16/B/D7	ICU 16/B/D7	5150542	1,20	14543912	0,93	1,29
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			3973905	0,92	13555890	0,87	

OPA1-S (lower band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 1/ set A

cal. mean

5465250

cal. mean

41570172

	calibrator		5782680	1,06	42065400	1,01	1,05
A1	H1/B	Control 17/B	6501600	1,19	42967628	1,03	1,15
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	36839820	0,89	#####
A4	Sample X/T	Sample X/T	excl.	#####	32757912	0,79	#####
A5	H3/B	Control 3/B	6432660	1,18	35733870	0,86	1,37
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	7877835	1,44	28368765	0,68	2,11
A8	CR4/B/D7	ICU 4/B/D7	6490800	1,19	30026880	0,72	1,64
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	5926095	1,08	26405840	0,64	1,71
A12	CR8/B/D7	ICU 11/B/D7	5904135	1,08	28186268	0,68	1,59
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	5565834	1,02	32125544	0,77	1,32
A16	CR10/B/D7	ICU 13/B/D7	7017840	1,28	37040626	0,89	1,44
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	8033212	1,47	34330306	0,83	1,78
A20	CR11/B/D7	ICU 9/B/D7	5454206	1,00	33386188	0,80	1,24
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		5147820	0,94	41074944	0,99	0,95

OPA1-S (lower band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean

6661993

cal. mean

39443078,5

	calibrator		5960070	0,89	31672407	0,80	1,11
B1	H4/B	Control 12/B	6200325	0,93	34761556	0,88	1,06
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	6801525	1,02	31297500	0,79	1,29
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	7059060	1,06	26830935	0,68	1,56
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	25290990	0,64	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	23090160	0,59	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25183170	0,64	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22048708	0,56	#####
B11	CR9/B/D0	ICU 3/B/D0	5120100	0,77	28875120	0,73	1,05
B12	CR9/B/D7	ICU 3/B/D7	4548150	0,68	32768164	0,83	0,82
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	5397747	0,81	32436270	0,82	0,99
B16	CR12/B/D7	ICU 20/B/D7	7798516	1,17	33377085	0,85	1,38
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	9759992	1,47	35632337	0,90	1,62
B20	CR13/B/D7	ICU 22/B/D7	11361976	1,71	35268420	0,89	1,91
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		7363916	1,11	47213750	1,20	0,92

OPA1-S (lower band)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 3/ set C		cal. mean	1345471,5	cal. mean	22029750		
	calibrator		1737731	1,29	20662800	0,94	1,38
C1	H15/B	Control 2/B	1589116	1,18	20219460	0,92	1,29
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	1797268	1,34	15286365	0,69	1,93
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	1386840	1,03	14141864	0,64	1,61
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	1674998	1,24	14270928	0,65	1,92
C8	CR3/B/D7	ICU 5/B/D7	1804488	1,34	13239864	0,60	2,23
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	1040474	0,77	18034032	0,82	0,94
C12	CR5/B/D7	ICU 8/B/D7	1137810	0,85	15358512	0,70	1,21
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	1766160	1,31	17619684	0,80	1,64
C16	CR7/B/D7	ICU 6/B/D7	1928182	1,43	17926524	0,81	1,76
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	1964025	1,46	14741100	0,67	2,18
C20	CR20/B/D7	ICU 7/B/D7	1964568	1,46	15625234	0,71	2,06
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		953212	0,71	23396700	1,06	0,67

OPA1-S (lower band)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	1052190		cal. mean	27123350	
calibrator			1299420	1,23	21501050	0,79	1,56
D1	H2/B	Control 4/B	1446660	1,37	15463789	0,57	2,41
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	467865	0,44	15370470	0,57	0,78
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	749385	0,71	13315368	0,49	1,45
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	747405	0,71	15176872	0,56	1,27
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	1577970	1,50	12146022	0,45	3,35
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	1109340	1,05	13352832	0,49	2,14
D12	CR2/B/D7	ICU 14/B/D7	1198350	1,14	14137728	0,52	2,19
D13	CR2/B/D180	ICU 14/B/D180	1285380	1,22	11800978	0,44	2,81
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	2140785	2,03	13229508	0,49	4,17
D18	CR17/B/D7	ICU 23/B/D7	1423890	1,35	11261786	0,42	3,26
D19	CR17/B/D180	ICU 23/B/D180	834750	0,79	10989998	0,41	1,96
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			804960	0,77	32745650	1,21	0,63

OPA1-S (lower band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 5/ set E
cal. mean
1124992
cal. mean
13612646,5

	calibrator		1314167	1,17	13056250	0,96	1,22
E1	H12/B	Control 7/B	2876588	2,56	15730101	1,16	2,21
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	890284	0,79	17836110	1,31	0,60
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	3802676	3,38	17137798	1,26	2,68
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	2508155	2,23	16121115	1,18	1,88
E8	CR6/B/D7	ICU 21/B/D7	2714850	2,41	12428020	0,91	2,64
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	2091915	1,86	12644640	0,93	2,00
E12	CR14/B/D7	ICU 10/B/D7	2358648	2,10	13423200	0,99	2,13
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####
E19	CR19/B/D0	ICU 15/B/D0	1357219	1,21	18376065	1,35	0,89
E20	CR19/B/D7	ICU 15/B/D7	2092420	1,86	18750060	1,38	1,35
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		935817	0,83	14169043	1,04	0,80

OPA1-S (lower band)

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 6/ set F
cal. mean
1426027
cal. mean
24950732

	calibrator		1944296	1,36	25509639	1,02	1,33
F1	H9/B	Control 9/B	2221267	1,56	32156010	1,29	1,21
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	2624670	1,84	31970444	1,28	1,44
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	1387957	0,97	25759862	1,03	0,94
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	1557990	1,09	24703341	0,99	1,10
F8	CR1/B/D7	ICU 1/B/D7	2589570	1,82	22559436	0,90	2,01
F9	CR1/B/D180	ICU 1/B/D180	2190382	1,54	26013983	1,04	1,47
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	2827440	1,98	24749871	0,99	2,00
F14	CR18/B/D7	ICU 18/B/D7	3493275	2,45	21148778	0,85	2,89
F15	CR18/B/D180	ICU 18/B/D180	2266110	1,59	21108311	0,85	1,88
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	1371460	0,96	27652779	1,11	0,87
F20	CR15/B/D7	ICU 12/B/D7	1403391	0,98	28943200	1,16	0,85
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		907758	0,64	24391825	0,98	0,65

OPA1-S (lower band)

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	2279205	cal. mean	15623570		
calibrator			2556360	1,12	17691250	1,13	
E15	CR16/B/D0	ICU 16/B/D0	3539884	1,55	14143732	0,91	1,72
E16	CR16/B/D7	ICU 16/B/D7	3578345	1,57	14543912	0,93	1,69
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			2002050	0,88	13555890	0,87	

OPA1-L vs. OPA1-S

WB ID		Control/Patient ID	OPA1-L density	OPA1-L normalized to calibrator	OPA1-S density	OPA1-S normalized to calibrator	OPA1-L/OPA1-S
GEL 1/ set A		cal. mean	6214612,5		cal. mean	5465250	
	calibrator		7497630	1,21	5782680	1,06	1,14
A1	H1/B	Control 17/B	7544340	1,21	6501600	1,19	1,02
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	5757120	1,05	#####
A4	Sample X/T	Sample X/T	excl.	#####	3401955	0,62	#####
A5	H3/B	Control 3/B	7926660	1,28	6432660	1,18	1,08
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	5397570	0,87	7877835	1,44	0,60
A8	CR4/B/D7	ICU 4/B/D7	5267475	0,85	6490800	1,19	0,71
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	5111010	0,82	5926095	1,08	0,76
A12	CR8/B/D7	ICU 11/B/D7	5631615	0,91	5904135	1,08	0,84
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	5443843	0,88	5565834	1,02	0,86
A16	CR10/B/D7	ICU 13/B/D7	5619645	0,90	7017840	1,28	0,70
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	5215408	0,84	8033212	1,47	0,57
A20	CR11/B/D7	ICU 9/B/D7	6976578	1,12	5454206	1,00	1,12
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		4931595	0,79	5147820	0,94	0,84

OPA1-L vs. OPA1-S

WB ID	Control/Patient ID	OPA1-L density	OPA1-L normalized to calibrator	OPA1-S density	OPA1-S normalized to calibrator	OPA1-L/OPA1-S
-------	--------------------	-------------------	---------------------------------------	-------------------	---------------------------------------	---------------

GEL 2/ set B

cal. mean

6667057

cal. mean

6661993

	calibrator		6068835	0,91	5960070	0,89	1,02
B1	H4/B	Control 12/B	5838840	0,88	6200325	0,93	0,94
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	5453280	0,82	6801525	1,02	0,80
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	5134680	0,77	7059060	1,06	0,73
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	5277105	0,79	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	4347990	0,65	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	5321070	0,80	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	4737735	0,71	#####
B11	CR9/B/D0	ICU 3/B/D0	3499245	0,52	5120100	0,77	0,68
B12	CR9/B/D7	ICU 3/B/D7	3088710	0,46	4548150	0,68	0,68
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	3439484	0,52	5397747	0,81	0,64
B16	CR12/B/D7	ICU 20/B/D7	4561744	0,68	7798516	1,17	0,58
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	5225616	0,78	9759992	1,47	0,54
B20	CR13/B/D7	ICU 22/B/D7	5349544	0,80	11361976	1,71	0,47
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		7265279	1,09	7363916	1,11	0,99

OPA1-L vs. OPA1-S

WB ID		Control/Patient ID	OPA1-L density	OPA1-L normalized to calibrator	OPA1-S density	OPA1-S normalized to calibrator	OPA1-L/OPA1-S
GEL 3/ set C		cal. mean	12017075		cal. mean	1345471,5	
	calibrator		11766732	0,98	1737731	1,29	0,76
C1	H15/B	Control 2/B	11093084	0,92	1589116	1,18	0,78
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	11870276	0,99	1797268	1,34	0,74
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	11324460	0,94	1386840	1,03	0,91
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	7537652	0,63	1674998	1,24	0,50
C8	CR3/B/D7	ICU 5/B/D7	9888942	0,82	1804488	1,34	0,61
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	12336970	1,03	1040474	0,77	1,33
C12	CR5/B/D7	ICU 8/B/D7	10709352	0,89	1137810	0,85	1,05
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	12870968	1,07	1766160	1,31	0,82
C16	CR7/B/D7	ICU 6/B/D7	10077036	0,84	1928182	1,43	0,59
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	12286003	1,02	1964025	1,46	0,70
C20	CR20/B/D7	ICU 7/B/D7	14452924	1,20	1964568	1,46	0,82
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		12267418	1,02	953212	0,71	1,44

OPA1-L vs. OPA1-S

WB ID	Control/Patient ID	OPA1-L density	OPA1-L normalized to calibrator	OPA1-S density	OPA1-S normalized to calibrator	OPA1-L/OPA1-S
-------	--------------------	-------------------	---------------------------------------	-------------------	---------------------------------------	---------------

GEL 4/ set D
cal. mean
9487102,5
cal. mean
1052190

	calibrator		10049400	1,06	1299420	1,23	0,86
D1	H2/B	Control 4/B	12756555	1,34	1446660	1,37	0,98
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	5276070	0,56	467865	0,44	1,25
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	4978980	0,52	749385	0,71	0,74
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	5876325	0,62	747405	0,71	0,87
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	7436700	0,78	1577970	1,50	0,52
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	8344395	0,88	1109340	1,05	0,83
D12	CR2/B/D7	ICU 14/B/D7	7844400	0,83	1198350	1,14	0,73
D13	CR2/B/D180	ICU 14/B/D180	8646750	0,91	1285380	1,22	0,75
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	13831470	1,46	2140785	2,03	0,72
D18	CR17/B/D7	ICU 23/B/D7	11146500	1,17	1423890	1,35	0,87
D19	CR17/B/D180	ICU 23/B/D180	13643550	1,44	834750	0,79	1,81
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		8924805	0,94	804960	0,77	1,23

OPA1-L vs. OPA1-S

WB ID		Control/Patient ID	OPA1-L density	OPA1-L normalized to calibrator	OPA1-S density	OPA1-S normalized to calibrator	OPA1-L/OPA1-S
GEL 5/ set E		cal. mean	6622605,5		cal. mean	1124992	
calibrator			6786894	1,02	1314167	1,17	0,88
E1	H12/B	Control 7/B	9449632	1,43	2876588	2,56	0,56
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	6779204	1,02	890284	0,79	1,29
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	7881524	1,19	3802676	3,38	0,35
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	9265439	1,40	2508155	2,23	0,63
E8	CR6/B/D7	ICU 21/B/D7	5076990	0,77	2714850	2,41	0,32
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	6513435	0,98	2091915	1,86	0,53
E12	CR14/B/D7	ICU 10/B/D7	7567282	1,14	2358648	2,10	0,55
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	1987965	1,77	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	2193349	1,95	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	28341	0,03	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	23265	0,02	#####
E19	CR19/B/D0	ICU 15/B/D0	8358292	1,26	1357219	1,21	1,05
E20	CR19/B/D7	ICU 15/B/D7	7947588	1,20	2092420	1,86	0,65
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			6458317	0,98	935817	0,83	1,17

OPA1-L vs. OPA1-S

WB ID		Control/Patient ID	OPA1-L density	OPA1-L normalized to calibrator	OPA1-S density	OPA1-S normalized to calibrator	OPA1-L/OPA1-S
GEL 6/ set F		cal. mean	9168219,5		cal. mean	1426027	
	calibrator		9497525	1,04	1944296	1,36	0,76
F1	H9/B	Control 9/B	12429009	1,36	2221267	1,56	0,87
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	10357650	1,13	2624670	1,84	0,61
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	8423622	0,92	1387957	0,97	0,94
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	7570035	0,83	1557990	1,09	0,76
F8	CR1/B/D7	ICU 1/B/D7	5860538	0,64	2589570	1,82	0,35
F9	CR1/B/D180	ICU 1/B/D180	6842224	0,75	2190382	1,54	0,49
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	9948240	1,09	2827440	1,98	0,55
F14	CR18/B/D7	ICU 18/B/D7	10238950	1,12	3493275	2,45	0,46
F15	CR18/B/D180	ICU 18/B/D180	9353295	1,02	2266110	1,59	0,64
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	14508383	1,58	1371460	0,96	1,65
F20	CR15/B/D7	ICU 12/B/D7	11674457	1,27	1403391	0,98	1,29
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		8838914	0,96	907758	0,64	1,51

OPA1-L vs. OPA1-S

WB ID		Control/Patient ID	OPA1-L density	OPA1-L normalized to calibrator	OPA1-S density	OPA1-S normalized to calibrator	OPA1-L/OPA1-S
GEL 7/ set ICU 16		cal. mean	4300663,5	cal. mean	2279205		
calibrator			4627422	1,08	2556360	1,12	
E15	CR16/B/D0	ICU 16/B/D0	4714126	1,10	3539884	1,55	0,71
E16	CR16/B/D7	ICU 16/B/D7	5150542	1,20	3578345	1,57	0,76
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			3973905	0,92	2002050	0,88	

Lipin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	5621254	cal. mean	25207886		
calibrator			6824688	1,21	26237938	1,04	1,17
A1	H1/B	Control 17/B	6087080	1,08	25525559	1,01	1,07
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	20146728	0,80	#####
A4	Sample X/T	Sample X/T	excl.	#####	19557452	0,78	#####
A5	H3/B	Control 3/B	3010946	0,54	20382535	0,81	0,66
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	4707270	0,84	20062467	0,80	1,05
A8	CR4/B/D7	ICU 4/B/D7	5608748	1,00	19913007	0,79	1,26
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	7493395	1,33	16150422	0,64	2,08
A12	CR8/B/D7	ICU 11/B/D7	6137261	1,09	17634541	0,70	1,56
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	1748165	0,31	17586742	0,70	0,45
A16	CR10/B/D7	ICU 13/B/D7	3793331	0,67	18543098	0,74	0,92
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	1579218	0,28	20294868	0,81	0,35
A20	CR11/B/D7	ICU 9/B/D7	3763317	0,67	21152021	0,84	0,80
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			4417820	0,79	24177834	0,96	0,82

Lipin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 2/ set B		cal. mean	5028018		cal. mean	28899095,5	
calibrator			4658160	0,93	27854350	0,96	0,96
B1	H4/B	Control 12/B	3105546	0,62	28818080	1,00	0,62
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	2547277	0,51	29889932	1,03	0,49
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	6616238	1,32	23959672	0,83	1,59
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	29444064	1,02	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	27218925	0,94	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25503469	0,88	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22796640	0,79	#####
B11	CR9/B/D0	ICU 3/B/D0	6701326	1,33	26684496	0,92	1,44
B12	CR9/B/D7	ICU 3/B/D7	4381743	0,87	27844674	0,96	0,90
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	6431596	1,28	29392896	1,02	1,26
B16	CR12/B/D7	ICU 20/B/D7	4443018	0,88	29492584	1,02	0,87
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	1716861	0,34	28194210	0,98	0,35
B20	CR13/B/D7	ICU 22/B/D7	3280513	0,65	25151087	0,87	0,75
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			5397876	1,07	29943841	1,04	1,04

Lipin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 3/ set C		cal. mean	6716594,5	cal. mean	34079936		
calibrator			7445284	1,11	34138962	1,00	1,11
C1	H15/B	Control 2/B	6939606	1,03	29730870	0,87	1,18
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	3478198	0,52	30098518	0,88	0,59
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	6764300	1,01	20994664	0,62	1,63
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	5170455	0,77	22660589	0,66	1,16
C8	CR3/B/D7	ICU 5/B/D7	5875144	0,87	18025875	0,53	1,65
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	3571394	0,53	19994035	0,59	0,91
C12	CR5/B/D7	ICU 8/B/D7	10724532	1,60	17391984	0,51	3,13
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	3761696	0,56	22085112	0,65	0,86
C16	CR7/B/D7	ICU 6/B/D7	6569055	0,98	18388750	0,54	1,81
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	1398952	0,21	24852613	0,73	0,29
C20	CR20/B/D7	ICU 7/B/D7	2437908	0,36	25778115	0,76	0,48
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			5987905	0,89	34020910	1,00	0,89

Lipin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	5003392,5	cal. mean	35190284		
calibrator			4916385	0,98	37167787	1,06	0,93
D1	H2/B	Control 4/B	780999	0,16	30949908	0,88	0,18
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	2343264	0,47	31992571	0,91	0,52
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	3833106	0,77	24406066	0,69	1,10
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	5075376	1,01	22372376	0,64	1,60
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	474642	0,09	18229138	0,52	0,18
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	4864735	0,97	17679473	0,50	1,94
D12	CR2/B/D7	ICU 14/B/D7	3194955	0,64	18461459	0,52	1,22
D13	CR2/B/D180	ICU 14/B/D180	2140380	0,43	19948365	0,57	0,75
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	2402595	0,48	22814050	0,65	0,74
D18	CR17/B/D7	ICU 23/B/D7	1136898	0,23	24909295	0,71	0,32
D19	CR17/B/D180	ICU 23/B/D180	6528764	1,30	26712850	0,76	1,72
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			5090400	1,02	33212781	0,94	1,08

Lipin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	6462175		cal. mean	13612646,5	
calibrator			4820140	0,75	13056250	0,96	0,78
E1	H12/B	Control 7/B	2479590	0,38	15730101	1,16	0,33
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	2967750	0,46	17836110	1,31	0,35
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	9304875	1,44	17137798	1,26	1,14
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	3829905	0,59	16121115	1,18	0,50
E8	CR6/B/D7	ICU 21/B/D7	5019168	0,78	12428020	0,91	0,85
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	3202554	0,50	12644640	0,93	0,53
E12	CR14/B/D7	ICU 10/B/D7	7848360	1,21	13423200	0,99	1,23
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####
E19	CR19/B/D0	ICU 15/B/D0	1985421	0,31	18376065	1,35	0,23
E20	CR19/B/D7	ICU 15/B/D7	2376000	0,37	18750060	1,38	0,27
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			8104210	1,25	14169043	1,04	1,20

Lipin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	7388423,5		cal. mean	24950732	
calibrator			8051805	1,09	25509639	1,02	1,07
F1	H9/B	Control 9/B	1032562	0,14	32156010	1,29	0,11
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	9614832	1,30	31970444	1,28	1,02
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	6314670	0,85	25759862	1,03	0,83
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	4756444	0,64	24703341	0,99	0,65
F8	CR1/B/D7	ICU 1/B/D7	7245348	0,98	22559436	0,90	1,08
F9	CR1/B/D180	ICU 1/B/D180	6471806	0,88	26013983	1,04	0,84
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	4864244	0,66	24749871	0,99	0,66
F14	CR18/B/D7	ICU 18/B/D7	983928	0,13	21148778	0,85	0,16
F15	CR18/B/D180	ICU 18/B/D180	6265429	0,85	21108311	0,85	1,00
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	4003164	0,54	27652779	1,11	0,49
F20	CR15/B/D7	ICU 12/B/D7	5740672	0,78	28943200	1,16	0,67
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			6725042	0,91	24391825	0,98	0,93

Lipin 1

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	17547382,5	cal. mean	40515874		
calibrator			19029000	1,08	43007588	1,06	
E15	CR16/B/D0	ICU 16/B/D0	6545527	0,37	37231612	0,92	0,41
E16	CR16/B/D7	ICU 16/B/D7	6247381	0,36	33408164	0,82	0,43
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			16065765	0,92	38024160	0,94	

Fatty acid synthase

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	7572334	cal. mean	26252326		
calibrator			7854088	1,04	27774192	1,06	0,98
A1	H1/B	Control 17/B	4867676	0,64	26549130	1,01	0,64
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	23658300	0,90	#####
A4	Sample X/T	Sample X/T	excl.	#####	20212632	0,77	#####
A5	H3/B	Control 3/B	4953916	0,65	22203368	0,85	0,77
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	5220204	0,69	22318120	0,85	0,81
A8	CR4/B/D7	ICU 4/B/D7	5002800	0,66	21786380	0,83	0,80
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	5968644	0,79	24650043	0,94	0,84
A12	CR8/B/D7	ICU 11/B/D7	6194760	0,82	26404414	1,01	0,81
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	4507140	0,60	21712328	0,83	0,72
A16	CR10/B/D7	ICU 13/B/D7	5211888	0,69	23581122	0,90	0,77
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	4378352	0,58	23880825	0,91	0,64
A20	CR11/B/D7	ICU 9/B/D7	4767092	0,63	21025925	0,80	0,79
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			7290580	0,96	24730460	0,94	1,02

Fatty acid synthase

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 2/ set B		cal. mean	7240354	cal. mean	28899095,5		
calibrator			7942452	1,10	27854350	0,96	1,14
B1	H4/B	Control 12/B	6050057	0,84	28818080	1,00	0,84
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	5421606	0,75	29889932	1,03	0,72
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	4475150	0,62	23959672	0,83	0,75
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	29444064	1,02	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	27218925	0,94	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25503469	0,88	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22796640	0,79	#####
B11	CR9/B/D0	ICU 3/B/D0	3355806	0,46	26684496	0,92	0,50
B12	CR9/B/D7	ICU 3/B/D7	2887364	0,40	27844674	0,96	0,41
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	4847777	0,67	29392896	1,02	0,66
B16	CR12/B/D7	ICU 20/B/D7	5014828	0,69	29492584	1,02	0,68
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	4528734	0,63	28194210	0,98	0,64
B20	CR13/B/D7	ICU 22/B/D7	5336940	0,74	25151087	0,87	0,85
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			6538256	0,90	29943841	1,04	0,87

Fatty acid synthase

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 3/ set C		cal. mean	8350909	cal. mean	22029750	
calibrator		8439712	1,01	20662800	0,94	1,08
C1 H15/B	Control 2/B	6847928	0,82	20219460	0,92	0,89
C2 H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3 H16/B	Control 14/B	3803608	0,46	15286365	0,69	0,66
C4 H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5 H17/B	Control 10/B	3776876	0,45	14141864	0,64	0,70
C6 H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7 CR3/B/D0	ICU 5/B/D0	3882956	0,46	14270928	0,65	0,72
C8 CR3/B/D7	ICU 5/B/D7	3442795	0,41	13239864	0,60	0,69
C9 CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10 CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11 CR5/B/D0	ICU 8/B/D0	2761804	0,33	18034032	0,82	0,40
C12 CR5/B/D7	ICU 8/B/D7	3843708	0,46	15358512	0,70	0,66
C13 CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14 CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15 CR7/B/D0	ICU 6/B/D0	3891830	0,47	17619684	0,80	0,58
C16 CR7/B/D7	ICU 6/B/D7	4403616	0,53	17926524	0,81	0,65
C17 CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18 CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19 CR20/B/D0	ICU 7/B/D0	3199284	0,38	14741100	0,67	0,57
C20 CR20/B/D7	ICU 7/B/D7	4282738	0,51	15625234	0,71	0,72
C21 CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22 CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator		8262106	0,99	23396700	1,06	0,93

Fatty acid synthase

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	16248662,5	cal. mean	27123350		
calibrator			14190797	0,87	21501050	0,79	1,10
D1	H2/B	Control 4/B	5799285	0,36	15463789	0,57	0,63
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	7191062	0,44	15370470	0,57	0,78
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	5841748	0,36	13315368	0,49	0,73
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	7754104	0,48	15176872	0,56	0,85
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	5544585	0,34	12146022	0,45	0,76
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	6563434	0,40	13352832	0,49	0,82
D12	CR2/B/D7	ICU 14/B/D7	5413804	0,33	14137728	0,52	0,64
D13	CR2/B/D180	ICU 14/B/D180	6278781	0,39	11800978	0,44	0,89
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	7807768	0,48	13229508	0,49	0,99
D18	CR17/B/D7	ICU 23/B/D7	7394760	0,46	11261786	0,42	1,10
D19	CR17/B/D180	ICU 23/B/D180	7196112	0,44	10989998	0,41	1,09
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			18306528	1,13	32745650	1,21	0,93

Fatty acid synthase

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH	
GEL 5/ set E		cal. mean	16457149,5		cal. mean	13612646,5		
calibrator			17163612	1,04	13056250	0,96	1,09	
E1	H12/B	Control 7/B	d.f.	#####	15730101	1,16	#####	
E2	H12/T	Control 7/T	d.f.	#####	d.n.r.	d.n.r.	#####	
E3	H13/B	Control 1/B	d.f.	#####	17836110	1,31	#####	
E4	H13/T	Control 1/T	d.f.	#####	d.n.r.	d.n.r.	#####	
E5	H18/B	Control 11/B		8570184	0,52	17137798	1,26	0,41
E6	H18/T	Control 11/T		d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0		9086238	0,55	16121115	1,18	0,47
E8	CR6/B/D7	ICU 21/B/D7		9156258	0,56	12428020	0,91	0,61
E9	CR6/T/D0	ICU 21/T/D0		d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7		d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0		7751099	0,47	12644640	0,93	0,51
E12	CR14/B/D7	ICU 10/B/D7		7514595	0,46	13423200	0,99	0,46
E13	CR14/T/D0	ICU 10/T/D0		d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7		d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####	
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####	
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####	
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####	
E19	CR19/B/D0	ICU 15/B/D0	d.f.	#####	18376065	1,35	#####	
E20	CR19/B/D7	ICU 15/B/D7	d.f.	#####	18750060	1,38	#####	
E21	CR19/T/D0	ICU 15/T/D0	d.f.	#####	d.n.r.	d.n.r.	#####	
E22	CR19/T/D7	ICU 15/T/D7	d.f.	#####	d.n.r.	d.n.r.	#####	
calibrator			15750687	0,96	14169043	1,04	0,92	

Fatty acid synthase

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	13593580	cal. mean	24950732		
calibrator			12678424	0,93	25509639	1,02	0,91
F1	H9/B	Control 9/B	8863360	0,65	32156010	1,29	0,51
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	11064284	0,81	31970444	1,28	0,64
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	10562640	0,78	25759862	1,03	0,75
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	6343304	0,47	24703341	0,99	0,47
F8	CR1/B/D7	ICU 1/B/D7	8683752	0,64	22559436	0,90	0,71
F9	CR1/B/D180	ICU 1/B/D180	9093568	0,67	26013983	1,04	0,64
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	6780928	0,50	24749871	0,99	0,50
F14	CR18/B/D7	ICU 18/B/D7	7652788	0,56	21148778	0,85	0,66
F15	CR18/B/D180	ICU 18/B/D180	7789936	0,57	21108311	0,85	0,68
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	6177600	0,45	27652779	1,11	0,41
F20	CR15/B/D7	ICU 12/B/D7	9864272	0,73	28943200	1,16	0,63
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			14508736	1,07	24391825	0,98	1,09

Fatty acid synthase

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	6403522	cal. mean	40515874		
calibrator			5669224	0,89	43007588	1,06	
E15	CR16/B/D0	ICU 16/B/D0	6573308	1,03	37231612	0,92	1,12
E16	CR16/B/D7	ICU 16/B/D7	7208890	1,13	33408164	0,82	1,37
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			7137820	1,11	38024160	0,94	

P-ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	9713542,5	cal. mean	26252326		
calibrator			9109665	0,94	27774192	1,06	0,89
A1	H1/B	Control 17/B	7832610	0,81	26549130	1,01	0,80
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	23658300	0,90	#####
A4	Sample X/T	Sample X/T	excl.	#####	20212632	0,77	#####
A5	H3/B	Control 3/B	6025230	0,62	22203368	0,85	0,73
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	6011820	0,62	22318120	0,85	0,73
A8	CR4/B/D7	ICU 4/B/D7	5617170	0,58	21786380	0,83	0,70
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	8766900	0,90	24650043	0,94	0,96
A12	CR8/B/D7	ICU 11/B/D7	7408350	0,76	26404414	1,01	0,76
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	3655440	0,38	21712328	0,83	0,46
A16	CR10/B/D7	ICU 13/B/D7	4927095	0,51	23581122	0,90	0,56
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	4401924	0,45	23880825	0,91	0,50
A20	CR11/B/D7	ICU 9/B/D7	3338415	0,34	21025925	0,80	0,43
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			10317420	1,06	24730460	0,94	1,13

P-ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 2/ set B		cal. mean	7870924	cal. mean	28899095,5		
calibrator			9334136	1,19	27854350	0,96	1,23
B1	H4/B	Control 12/B	9167914	1,16	28818080	1,00	1,17
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	6240768	0,79	29889932	1,03	0,77
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	6113138	0,78	23959672	0,83	0,94
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	29444064	1,02	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	27218925	0,94	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	25503469	0,88	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	22796640	0,79	#####
B11	CR9/B/D0	ICU 3/B/D0	5501979	0,70	26684496	0,92	0,76
B12	CR9/B/D7	ICU 3/B/D7	7260090	0,92	27844674	0,96	0,96
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	7864785	1,00	29392896	1,02	0,98
B16	CR12/B/D7	ICU 20/B/D7	8378937	1,06	29492584	1,02	1,04
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	6854256	0,87	28194210	0,98	0,89
B20	CR13/B/D7	ICU 22/B/D7	5650169	0,72	25151087	0,87	0,82
B21	CR13/T/D0	ICU 22/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B22	CR13/T/D7	ICU 22/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			6407712	0,81	29943841	1,04	0,79

P-ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 3/ set C		cal. mean	18353954,5	cal. mean	22029750		
calibrator			15920149	0,87	20662800	0,94	0,92
C1	H15/B	Control 2/B	10633891	0,58	20219460	0,92	0,63
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	4900194	0,27	15286365	0,69	0,38
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	7010075	0,38	14141864	0,64	0,59
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	5969389	0,33	14270928	0,65	0,50
C8	CR3/B/D7	ICU 5/B/D7	6788840	0,37	13239864	0,60	0,62
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	6184819	0,34	18034032	0,82	0,41
C12	CR5/B/D7	ICU 8/B/D7	11138389	0,61	15358512	0,70	0,87
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	4595216	0,25	17619684	0,80	0,31
C16	CR7/B/D7	ICU 6/B/D7	6394745	0,35	17926524	0,81	0,43
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	2881302	0,16	14741100	0,67	0,23
C20	CR20/B/D7	ICU 7/B/D7	6794989	0,37	15625234	0,71	0,52
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			20787760	1,13	23396700	1,06	1,07

P-ACLY

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 4/ set D

cal. mean

12835182,5

cal. mean

27123350

	calibrator		14358274	1,12	21501050	0,79	1,41
D1	H2/B	Control 4/B	2416834	0,19	15463789	0,57	0,33
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	7152448	0,56	15370470	0,57	0,98
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	5201685	0,41	13315368	0,49	0,83
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	6834082	0,53	15176872	0,56	0,95
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	3408186	0,27	12146022	0,45	0,59
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	4699244	0,37	13352832	0,49	0,74
D12	CR2/B/D7	ICU 14/B/D7	3572611	0,28	14137728	0,52	0,53
D13	CR2/B/D180	ICU 14/B/D180	3557224	0,28	11800978	0,44	0,64
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	3517536	0,27	13229508	0,49	0,56
D18	CR17/B/D7	ICU 23/B/D7	3981198	0,31	11261786	0,42	0,75
D19	CR17/B/D180	ICU 23/B/D180	5116810	0,40	10989998	0,41	0,98
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		11312091	0,88	32745650	1,21	0,73

P-ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	10704106	cal. mean	13612646,5		
calibrator			9021428	0,84	13056250	0,96	0,88
E1	H12/B	Control 7/B	1431728	0,13	15730101	1,16	0,12
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	1268036	0,12	17836110	1,31	0,09
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	6358304	0,59	17137798	1,26	0,47
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	5236502	0,49	16121115	1,18	0,41
E8	CR6/B/D7	ICU 21/B/D7	5252188	0,49	12428020	0,91	0,54
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	11309100	1,06	12644640	0,93	1,14
E12	CR14/B/D7	ICU 10/B/D7	13765684	1,29	13423200	0,99	1,30
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	Sample W/B/D0	Sample W/B/D0	excl.	#####	13286318	0,98	#####
E16	Sample W/B/D7	Sample W/B/D7	excl.	#####	14503083	1,07	#####
E17	Sample W/T/D0	Sample W/T/D0	excl.	#####	15358237	1,13	#####
E18	Sample W/T/D7	Sample W/T/D7	excl.	#####	14918082	1,10	#####
E19	CR19/B/D0	ICU 15/B/D0	5350444	0,50	18376065	1,35	0,37
E20	CR19/B/D7	ICU 15/B/D7	5886630	0,55	18750060	1,38	0,40
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			12386784	1,16	14169043	1,04	1,11

P-ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 6/ set F		cal. mean	16038867,5	cal. mean	24950732		
calibrator			14696712	0,92	25509639	1,02	0,90
F1	H9/B	Control 9/B	6472752	0,40	32156010	1,29	0,31
F2	H9/T	Control 9/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F3	H10/B	Control 15/B	9608445	0,60	31970444	1,28	0,47
F4	H10/T	Control 15/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F5	H14/B	Control 20/B	8116636	0,51	25759862	1,03	0,49
F6	H14/T	Control 20/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F7	CR1/B/D0	ICU 1/B/D0	2297435	0,14	24703341	0,99	0,14
F8	CR1/B/D7	ICU 1/B/D7	7470650	0,47	22559436	0,90	0,52
F9	CR1/B/D180	ICU 1/B/D180	8836517	0,55	26013983	1,04	0,53
F10	CR1/T/D0	ICU 1/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F11	CR1/T/D7	ICU 1/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F12	CR1/T/D180	ICU 1/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F13	CR18/B/D0	ICU 18/B/D0	4034045	0,25	24749871	0,99	0,25
F14	CR18/B/D7	ICU 18/B/D7	4231316	0,26	21148778	0,85	0,31
F15	CR18/B/D180	ICU 18/B/D180	12209190	0,76	21108311	0,85	0,90
F16	CR18/T/D0	ICU 18/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F17	CR18/T/D7	ICU 18/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F18	CR18/T/D180	ICU 18/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F19	CR15/B/D0	ICU 12/B/D0	9540248	0,59	27652779	1,11	0,54
F20	CR15/B/D7	ICU 12/B/D7	5780017	0,36	28943200	1,16	0,31
F21	CR15/T/D0	ICU 12/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
F22	CR15/T/D7	ICU 12/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			17381023	1,08	24391825	0,98	1,11

P-ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 7/ set ICU 16		cal. mean	5863873	cal. mean	40515874		
calibrator			5884182	1,00	43007588	1,06	
E15	CR16/B/D0	ICU 16/B/D0	4460256	0,76	37231612	0,92	0,83
E16	CR16/B/D7	ICU 16/B/D7	3760868	0,64	33408164	0,82	0,78
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			5843564	1,00	38024160	0,94	

ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 1/ set A		cal. mean	5194952	cal. mean	32397000		
calibrator			3387750	0,65	34486272	1,06	0,61
A1	H1/B	Control 17/B	4945293	0,95	36056208	1,11	0,86
A2	H1/T	Control 17/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A3	Sample X/B	Sample X/B	excl.	#####	27900540	0,86	#####
A4	Sample X/T	Sample X/T	excl.	#####	28724784	0,89	#####
A5	H3/B	Control 3/B	2473812	0,48	28793205	0,89	0,54
A6	H3/T	Control 3/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A7	CR4/B/D0	ICU 4/B/D0	5037407	0,97	24152530	0,75	1,30
A8	CR4/B/D7	ICU 4/B/D7	3137814	0,60	22376194	0,69	0,87
A9	CR4/T/D0	ICU 4/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A10	CR4/T/D7	ICU 4/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A11	CR8/B/D0	ICU 11/B/D0	9684073	1,86	23807098	0,73	2,54
A12	CR8/B/D7	ICU 11/B/D7	11239345	2,16	24407424	0,75	2,87
A13	CR8/T/D0	ICU 11/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A14	CR8/T/D7	ICU 11/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A15	CR10/B/D0	ICU 13/B/D0	2203245	0,42	23284992	0,72	0,59
A16	CR10/B/D7	ICU 13/B/D7	3645226	0,70	23631488	0,73	0,96
A17	CR10/T/D0	ICU 13/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A18	CR10/T/D7	ICU 13/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A19	CR11/B/D0	ICU 9/B/D0	1725120	0,33	33965550	1,05	0,32
A20	CR11/B/D7	ICU 9/B/D7	1783744	0,34	32601396	1,01	0,34
A21	CR11/T/D0	ICU 9/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
A22	CR11/T/D7	ICU 9/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			7002154	1,35	30307728	0,94	1,44

ACLY

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 2/ set B

cal. mean	5192216	cal. mean	27865678,5
-----------	---------	-----------	------------

	calibrator		5467104	1,05	30163392	1,08	0,97
B1	H4/B	Control 12/B	7437985	1,43	27105885	0,97	1,47
B2	H4/T	Control 12/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B3	H5/B	Control 6/B	2197800	0,42	28410290	1,02	0,42
B4	H5/T	Control 6/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B5	H6/B	Control 8/B	3086349	0,59	20371736	0,73	0,81
B6	H6/T	Control 8/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B7	Sample Y/B	Sample Y/B	excl.	#####	22090275	0,79	#####
B8	Sample Y/T	Sample Y/T	excl.	#####	19395114	0,70	#####
B9	Sample Z/B	Sample Z/B	excl.	#####	20887034	0,75	#####
B10	Sample Z/T	Sample Z/T	excl.	#####	18124046	0,65	#####
B11	CR9/B/D0	ICU 3/B/D0	3617246	0,70	21057888	0,76	0,92
B12	CR9/B/D7	ICU 3/B/D7	4833000	0,93	20873640	0,75	1,24
B13	CR9/T/D0	ICU 3/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B14	CR9/T/D7	ICU 3/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B15	CR12/B/D0	ICU 20/B/D0	4519845	0,87	26873664	0,96	0,90
B16	CR12/B/D7	ICU 20/B/D7	5448193	1,05	23672428	0,85	1,24
B17	CR12/T/D0	ICU 20/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B18	CR12/T/D7	ICU 20/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
B19	CR13/B/D0	ICU 22/B/D0	4276890	0,82	22632840	0,81	1,01
B20	CR13/B/D7	ICU 22/B/D7	d.f.	#####	20968515	0,75	#####
B21	CR13/T/D0	ICU 22/T/D0	d.f.	#####	d.n.r.	d.n.r.	#####
B22	CR13/T/D7	ICU 22/T/D7	d.f.	#####	d.n.r.	d.n.r.	#####
	calibrator		4917328	0,95	25567965	0,92	1,03

ACLY

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 3/ set C

cal. mean	5477333	cal. mean	33101325,5
-----------	---------	-----------	------------

	calibrator		5443916	0,99	39634150	1,20	0,83
C1	H15/B	Control 2/B	6554208	1,20	41027124	1,24	0,97
C2	H15/T	Control 2/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C3	H16/B	Control 14/B	3664314	0,67	36119890	1,09	0,61
C4	H16/T	Control 14/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C5	H17/B	Control 10/B	6015742	1,10	32155610	0,97	1,13
C6	H17/T	Control 10/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C7	CR3/B/D0	ICU 5/B/D0	2926272	0,53	31662000	0,96	0,56
C8	CR3/B/D7	ICU 5/B/D7	2895705	0,53	32580420	0,98	0,54
C9	CR3/T/D0	ICU 5/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C10	CR3/T/D7	ICU 5/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C11	CR5/B/D0	ICU 8/B/D0	1331088	0,24	35399622	1,07	0,23
C12	CR5/B/D7	ICU 8/B/D7	2021404	0,37	28886670	0,87	0,42
C13	CR5/T/D0	ICU 8/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C14	CR5/T/D7	ICU 8/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C15	CR7/B/D0	ICU 6/B/D0	209367	0,04	25648220	0,77	0,05
C16	CR7/B/D7	ICU 6/B/D7	516810	0,09	24801885	0,75	0,13
C17	CD7/T/D0	ICU 6/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C18	CR7/T/D7	ICU 6/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C19	CR20/B/D0	ICU 7/B/D0	468424	0,09	23268402	0,70	0,12
C20	CR20/B/D7	ICU 7/B/D7	1225900	0,22	26161994	0,79	0,28
C21	CR20/T/D0	ICU 7/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
C22	CR20/T/D7	ICU 7/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
	calibrator		5510750	1,01	26568501	0,80	1,25

ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 4/ set D		cal. mean	5435325		cal. mean	37613053,5	
calibrator			4352572	0,80	35832900	0,95	0,84
D1	H2/B	Control 4/B	69930	0,01	33954570	0,90	0,01
D2	H2/T	Control 4/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D3	H7/B	Control 13/B	4406524	0,81	37204616	0,99	0,82
D4	H7/T	Control 13/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D5	H8/B	Control 18/B	7407610	1,36	27370416	0,73	1,87
D6	H8/T	Control 18/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D7	H11/B	Control 19/B	6382086	1,17	33493056	0,89	1,32
D8	H11/T	Control 19/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D9	H19/B	Control 5/B	3019946	0,56	29291952	0,78	0,71
D10	H19/T	Control 5/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D11	CR2/B/D0	ICU 14/B/D0	4832530	0,89	26753646	0,71	1,25
D12	CR2/B/D7	ICU 14/B/D7	4020354	0,74	28843383	0,77	0,96
D13	CR2/B/D180	ICU 14/B/D180	3108674	0,57	26981280	0,72	0,80
D14	CR2/T/D0	ICU 14/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D15	CR2/T/D7	ICU 14/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D16	CR2/T/D180	ICU 14/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D17	CR17/B/D0	ICU 23/B/D0	7653903	1,41	28709933	0,76	1,84
D18	CR17/B/D7	ICU 23/B/D7	4378044	0,81	28537335	0,76	1,06
D19	CR17/B/D180	ICU 23/B/D180	5081850	0,93	29401155	0,78	1,20
D20	CR17/T/D0	ICU 23/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D21	CR17/T/D7	ICU 23/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
D22	CR17/T/D180	ICU 23/T/D180	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			6518078	1,20	39393207	1,05	1,15

ACLY

WB ID		Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
GEL 5/ set E		cal. mean	5039923		cal. mean	33459356	
calibrator			3320096	0,66	32367400	0,97	0,68
E1	H12/B	Control 7/B	1533465	0,30	33789330	1,01	0,30
E2	H12/T	Control 7/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E3	H13/B	Control 1/B	4064974	0,81	32255844	0,96	0,84
E4	H13/T	Control 1/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E5	H18/B	Control 11/B	9602429	1,91	27424188	0,82	2,32
E6	H18/T	Control 11/T	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E7	CR6/B/D0	ICU 21/B/D0	5665320	1,12	29010256	0,87	1,30
E8	CR6/B/D7	ICU 21/B/D7	4850013	0,96	22604584	0,68	1,42
E9	CR6/T/D0	ICU 21/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E10	CR6/T/D7	ICU 21/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E11	CR14/B/D0	ICU 10/B/D0	3270876	0,65	24246002	0,72	0,90
E12	CR14/B/D7	ICU 10/B/D7	4478220	0,89	26351913	0,79	1,13
E13	CR14/T/D0	ICU 10/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E14	CR14/T/D7	ICU 10/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E15	CR16/B/D0	ICU 16/B/D0	4027584	0,80	29617830	0,89	0,90
E16	CR16/B/D7	ICU 16/B/D7	3735952	0,74	26893768	0,80	0,92
E17	CR16/T/D0	ICU 16/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E18	CR16/T/D7	ICU 16/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E19	CR19/B/D0	ICU 15/B/D0	2686005	0,53	32525266	0,97	0,55
E20	CR19/B/D7	ICU 15/B/D7	3317625	0,66	38332904	1,15	0,57
E21	CR19/T/D0	ICU 15/T/D0	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
E22	CR19/T/D7	ICU 15/T/D7	d.n.r.	d.n.r.	d.n.r.	d.n.r.	d.n.r.
calibrator			6759750	1,34	34551312	1,03	1,30

ACLY

WB ID	Control/Patient ID	protein density	protein normalized to calibrator	GAPDH density	GAPDH normalized to calibrator	protein/GAPDH
-------	--------------------	-----------------	----------------------------------	---------------	--------------------------------	---------------

GEL 6/ set F

cal. mean	#####
-----------	-------

cal. mean	37954347
-----------	----------

calibrator			d.f.	#####	26550450	0,70	#####
F1	H9/B	Control 9/B	d.f.	#####	30236175	0,80	#####
F2	H9/T	Control 9/T	d.f.	#####	d.n.r.	d.n.r.	#####
F3	H10/B	Control 15/B	d.f.	#####	28345900	0,75	#####
F4	H10/T	Control 15/T	d.f.	#####	d.n.r.	d.n.r.	#####
F5	H14/B	Control 20/B	d.f.	#####	27331912	0,72	#####
F6	H14/T	Control 20/T	d.f.	#####	d.n.r.	d.n.r.	#####
F7	CR1/B/D0	ICU 1/B/D0	d.f.	#####	26053802	0,69	#####
F8	CR1/B/D7	ICU 1/B/D7	d.f.	#####	27329985	0,72	#####
F9	CR1/B/D180	ICU 1/B/D180	d.f.	#####	29696880	0,78	#####
F10	CR1/T/D0	ICU 1/T/D0	d.f.	#####	d.n.r.	d.n.r.	#####
F11	CR1/T/D7	ICU 1/T/D7	d.f.	#####	d.n.r.	d.n.r.	#####
F12	CR1/T/D180	ICU 1/T/D180	d.f.	#####	d.n.r.	d.n.r.	#####
F13	CR18/B/D0	ICU 18/B/D0	d.f.	#####	31147879	0,82	#####
F14	CR18/B/D7	ICU 18/B/D7	d.f.	#####	34315678	0,90	#####
F15	CR18/B/D180	ICU 18/B/D180	d.f.	#####	36800640	0,97	#####
F16	CR18/T/D0	ICU 18/T/D0	d.f.	#####	d.n.r.	d.n.r.	#####
F17	CR18/T/D7	ICU 18/T/D7	d.f.	#####	d.n.r.	d.n.r.	#####
F18	CR18/T/D180	ICU 18/T/D180	d.f.	#####	d.n.r.	d.n.r.	#####
F19	CR15/B/D0	ICU 12/B/D0	d.f.	#####	45890106	1,21	#####
F20	CR15/B/D7	ICU 12/B/D7	d.f.	#####	39184088	1,03	#####
F21	CR15/T/D0	ICU 12/T/D0	d.f.	#####	d.n.r.	d.n.r.	#####
F22	CR15/T/D7	ICU 12/T/D7	d.f.	#####	d.n.r.	d.n.r.	#####
calibrator				#####	49358244	1,30	#####