



mouse anti-ATXN2 (x1000),  
Hematoxylin



| Fixiation / preservation method |                                                     |                                                            |
|---------------------------------|-----------------------------------------------------|------------------------------------------------------------|
| 4% PFA-fixed, frozen floating   |                                                     | formalin-fixed, paraffin-embedded.                         |
| Heating preparation             |                                                     |                                                            |
| No preparation                  | Autoclaving (121°C, 20 min) in Citrate buffer, pH 6 | Pressure cooking (120°C, 10 min) in Tris-EDTA buffer, pH 9 |