

Supplemental Table S1: Acquisition Protocol

Site	Positronmed, Santiago, Chile	
PET/CT scanner	Biograph Vision 450 (6R/64CT), Siemens	Biograph mCT Flow 20-3R
CT reference (mAs)	190	170
CT peak kilovoltage (kV)	100	140
CT slice thickness (mm)	2	2
CT slice increment (mm)	1	2
PET reconstruction	OSEM algorithm (PSF & TOF)	OSEM algorithm (PSF & TOF)
Iterations	4	2
Subsets	5	21
Matrix	220 x 220	200 x 200
Corrections	Gaussian FWHM 4.0 mm	Gaussian FWHM 4.0 mm